For Research Use Only. Not for use in diagnostic procedures.
--

PURIFICATION that Maintains Protein Activity

c-Myc tagged Protein

MILD PURIFICATION GEL

(MoAb. clone 1G4)

CODE No. 3306, 3307

Purification with neutral pH (to maintain Activity and Conformation)

High Quality and Efficiency

Rapid, High-Yield Purification

Product Description

The ability to isolate and study a purified protein lies at the heart of modern biochemistry. Researchers in many fields require highly purified, active proteins for studies involving signaling pathways, enzymology, receptor binding, DNA binding, post-transcriptional modifications, and much more.

The method of purification is one of the important keys for maintaining protein structure and function. The c-Myc tagged Protein MILD PURIFICATION GEL is designed for the isolation of c-Myc tagged protein from cell lysates and culture supernatant. Severe conditions such as acidic or alkaline elution can denature protein structure while maintaining a neutral pH can preserve protein activity and conformation. MBL has developed the Anti-c-Myc tag Gel to quickly and efficiently purify c-Myc tagged proteins at neutral pH to maintain protein activity and native conformation. The c-Myc tag peptide competitively elutes the c-Myc tagged protein from Anti-c-Myc tag Gel while minimally interfering with protein function.

The affinity of the anti-c-Myc tag antibody is also very important. Very high affinity antibody cannot be dissociated from the c-Myc tagged protein, and low affinity antibody is not sufficient for binding. MBL has optimized the affinity of the anti-c-Myc tag antibody to permit high yields and efficient purification of functional c-Myc-tagged proteins for biochemical studies.

Anti-c-Myc tag monoclonal antibody (clone 1G4) used by this gel shows the best affinity to refine the c-Myc tag fusion protein under a mild condition.

Components

	Quantity	
	CODE No. 3306	CODE No. 3307
<u>Anti-c-Myc tag Gel</u>	1 mL × 1	1 mL × 5
50 % slurry: 1 mL Gel in 2 mL total volume in PBS with 0.1% ProClin 150 as preservative		
<u>Elution Peptide</u>	1 mg × 1	1 mg × 5
Lyophilized form: c-Myc tag peptide (EQKLISEEDL), reconstitute the Elution Peptide with 1 mL of distilled water. 1 mg in 1 mL PBS after reconstitution		

Product Capacity

The purification capacity of the Anti-c-Myc tag Gel varies depending upon the c-Myc tagged protein. For example, 1 mL of Anti-c-Myc Gel bound 1.0 mg of a c-Myc tagged protein (27 kDa) and eluted 0.64 mg of purified protein in our hands.

Storage

Store for up to 1 year from date of receipt at 2-8°C. Do not freeze.

Material Preparation

Prepare the following reagents before affinity purification.

Equilibration buffer	: PBS
Elution buffer	: 0.1-0.5 mg/mL c-myc tag peptide (EQKLISEEDL) in Equilibration buffer. Fully reconstitute the Elution peptide with 1 mL of distilled water, and mix it with 1-9 mL of Equilibration buffer in another tube. Store the reconstituted Elution peptide in aliquots at -20°C.
Regeneration buffer	: 0.17 M Glycine-HCl, (pH 2.3).
Column storage buffer	: PBS/Preservative (e.g. 0.1% ProClin 150 and 0.09% NaN ₃).

Protocols

- Notes:
1. Cellular debris and particulate matter must be removed prior to purification. The protein extract should be centrifuged ($10,000\text{-}20,000 \times g$ for 15 min) and filtered with a $0.45 \mu\text{m}$ filter to remove any remaining cells and particulates.
 2. Highly viscous samples containing chromosomal DNA or RNA should be sonicated or treated with nuclease to reduce viscosity.

A. Column preparation

1. Place the empty chromatography column on rack or stand.
2. Rinse the column with Equilibration buffer.
3. Resuspend Anti-c-Myc tag Gel by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
4. Transfer the desired volume to the column. Drain column storage buffer.
Do not allow the column to dry out.
5. Wash the resin with 10 bed-volumes of Regeneration buffer.
6. Immediately wash the resin with 10 bed-volumes of Equilibration buffer.

B. Loading column

1. Load the lysate on the top of the column under gravity flow.
Note: The binding efficiency depends on the fusion protein and on the loading flow rate.
The binding efficiency may be increased by multiple passes over the column.
2. Collect the flow-through into clean collection tubes
3. Wash the column with PBS until the OD₂₈₀ is <0.01

C. Elution


1. Elute the bound protein from the column with 5 bed-volumes of Elution buffer.
2. Collect the elution fraction into clean collection tubes.

D. Regeneration and storage

1. Wash the column with 10 bed-volumes of Regeneration buffer.
2. Immediately wash the column with 10 bed-volumes of Column storage buffer.
3. Store at $+2$ to $+8^{\circ}\text{C}$ in 2 bed-volumes of Column storage buffer.
Note: Poured columns containing the anti-c-Myc tag Gel may be used at least 10 times, depending on the usage conditions.

Example of Purification Results using c-Myc tagged Protein MILD PURIFICATION GEL

Elution under various conditions (using spin column.)


Lane 1: 0.1 mg/mL peptide
Lane 2: 1.0 mg/mL peptide
Lane 3: pH 4.0
Lane 4: pH 3.0
Lane 5: pH 2.7

The c-Myc tag fusion protein was purified under various conditions. The amount of eluted protein was nearly identical for all tested conditions excluding pH 4.0.

SDS-PAGE (Coomassie Brilliant Blue Staining)

Purification using open column



The c-Myc tagged protein was loaded onto the open column (bed vol. ~0.5 mL). Elution was carried out with 0.1 mg/mL peptide. Each fraction was 0.5 mL.

SDS-PAGE (Coomassie Brilliant Blue Staining)

Purification using HPLC system


N-terminus c-Myc-tagged GST

Column bed volume	: 1 mL
Sample	: c-Myc-tagged GST/293T (5×10^7 cells)
Cell lysis buffer	: 10 mM Tris-HCl, 150 mM NaCl, 1% NP-40 (pH 7.5)
Wash buffer	: 0.05% Tween/PBS (pH 7.2)
Elution buffer	: 0.5 mg/ml c-Myc peptide in PBS
Flow rate	: 0.75 mL/min
Chromatography system	: ÄKTA explorer 10S (GE Healthcare)

Purification of c-Myc tagged protein from mammalian cells


c-myc tagged SIRT1/Sir2 Deacetylase (N-terminal tagged, using spin column)

Human embryonic kidney cells (293T) were transfected and cultured for 70 hours. Cells were then lysed in the Lysis buffer (500 μ L/100 mm dish) and purified using a spin column. Elution was carried out with 0.1 mg/mL c-Myc tag peptide.

SDS-PAGE (Coomassie Brilliant Blue Staining)

c-Myc tagged GFRA3 (C-terminal tagged, using open column)

Human embryonic kidney cells (293T) were transfected and cultured for 70 hours. Culture supernatant was loaded onto the column (bed vol. 0.5 mL). Elution was carried out with 2.5 mL of 0.1 mg/mL c-Myc tag peptide. Each fraction was 0.5 mL.

SDS-PAGE (Coomassie Brilliant Blue Staining)

Additional Information

Several reagents were examined whether or not they were suitable for use with the c-Myc tagged Protein MILD PURIFICATION GEL. The results are listed below.

Chaotropic agents

Urea	1 M	Yes
Guanidine-HCl	1 M	No

Reducing agents

DTT	10 mM	Yes
2-Mercaptoethanol	10 mM	Yes

Surfactants

Nonionic	Tween-20	5%	Yes
	TritonX-100	5%	Yes
	NP40	5%	Yes
	Digitonin	1%	Yes
	n-Octyl-²-D-gulcoside	1%	Yes
Zwitterionic	CHAPS	1%	Yes
	CHAPSO	1%	Yes
Anionic	SDS	0.05%	No
	Sodium Deoxycholate	0.1%	Yes
	Sodium Deoxycholate	0.5%	No

Others

NaCl	1 M	Yes
Glycerol	10%	Yes
EDTA	10 mM	Yes

The “Yes” indicates the reagents can be used in the Lysis buffer for this Gel up to the indicated concentration. The “No” indicates the reagents cannot be used in the Lysis buffer for this Gel at the indicated concentration.

RELATED PRODUCTS:

3305	c-Myc-tagged Protein Mild Purification Kit ver.2	20 purifications
3305A	c-Myc-tagged Protein Mild Purification Kit ver.2	2 purifications
3340	c-Myc-tagged Protein Magnetic Purification Kit	1 kit
3340A	c-Myc-tagged Protein Magnetic Purification Kit (Trial Kit)	1 kit
3306	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 1 mg x 1	
3307	c-Myc-tagged Protein Mild Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 1 mg x 5	
3306K	c-Myc-tagged Protein PURIFICATION CARTRIDGE	1 mL x 1
3300-205	c-Myc-tag peptide	1 mg x 5
3310	His-tagged Protein Purification Kit	20 purifications
3310A	His-tagged Protein Purification Kit (Trial Kit)	2 purifications
3311	His-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 2 mg x 5	
3312	His-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 2 mg x 25	
3310-205	His-tag peptide	2 mg x 5
3317	V5-tagged Protein Purification Kit Ver.2	20 purifications
3317A	V5-tagged Protein Purification Kit Ver.2 (Trial Kit)	2 purifications
3341	V5-tagged Protein Magnetic Purification Kit	1 kit
3341A	V5-tagged Protein Magnetic Purification Kit (Trial Kit)	1 kit
3318	V5-tagged Protein Purification Gel Ver.2	Gel: 1 mL x 1
3315-205	V5-tag peptide	2 mg x 5
3320	HA-tagged Protein Purification Kit	20 purifications
3320A	HA-tagged Protein Purification Kit (Trial Kit)	2 purifications
3342	HA-tagged Protein Magnetic Purification Kit	1 kit
3342A	HA-tagged Protein Magnetic Purification Kit (Trial Kit)	1 kit
3321	HA-tagged Protein Purification Gel	1 mL
3320-205	HA-tag peptide	2 mg x 5
3325	DDDDK-tagged Protein Purification Kit	20 purifications
3325A	DDDDK-tagged Protein Purification Kit (Trial Kit)	2 purifications
3343	DDDDK-tagged Protein Magnetic Purification Kit	1 kit
3343A	DDDDK-tagged Protein Magnetic Purification Kit (Trial Kit)	1 kit
3326	DDDDK-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 1, Peptide: 1 mg x 5	
3327	DDDDK-tagged Protein Purification Gel with Elution Peptide	
	Gel: 1 mL x 5, Peptide: 1 mg x 25	
3328	DDDDK-tagged Protein Purification Gel	5 mL x 1
3329	DDDDK-tagged Protein Purification Gel	5 mL x 5
3326K	DDDDK-tagged Protein PURIFICATION CARTRIDGE	1 mL x 1
3325-205	DDDDK-tag peptide	1 mg x 5

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