For Research Use Only. Not for use in diagnostic procedures.

V5-tagged Protein PURIFICATION GEL Ver.2

(MoAb. clone OZA3)

Code No. 3318

*Purification to Maintain Protein Activity
from eukaryote cell lysate and culture supernatant*

Product Description

The ability to isolate and study a purified protein lies at the heart of modern biochemistry. Researchers in many fields require highly purified, active proteins for studies involving signaling pathways, enzymology, receptor binding, DNA binding, post-transcriptional modifications, and much more. Thus, choosing a method of purification is an important aspect in maintaining protein structure and function.

Recombinant tagged protein purification methods and kits are now widely recognized. The V5 epitope tag, GKPIPNPLLGLDST, is derived from the P/V proteins of paramyxovirus SV5. This V5-tag is also commonly used in such as mammalian and insect cell expression vectors.

MBL's V5-tagged Protein PURIFICATION GEL Ver.2 is designed for the isolation of V5-tagged protein from cell culture supernatants and cell lysate under neutral pH condition. Severe conditions such as acidic or alkaline elution denature protein structure. However, a neutral pH elution can preserve protein activity and native conformation. Therefore, MBL has developed the Anti-V5-tag Gel to purify V5-tagged proteins quickly and efficiently at neutral pH to maintain the activity and conformation of purified proteins. The elution of V5-tagged proteins from the Gel is achieved by the addition of the V5-tag peptide (GKPIPNPLLGLDST). As the V5-tag peptide competes with V5-tagged proteins on the Gel, the purified proteins do not lose the protein activity.

Components

	Quantity
	Code No. 3318
<u>Anti-V5-tag Gel</u>	1 mL × 1 vials

50% slurry: 1 mL Gel in 2 mL total volume in PBS with 0.1 % ProClin 150 as preservative.

Elution Peptide (GKPIPNPLLGLDST) is NOT included.

It is available as Code No. 3315-205 (2 mg × 5 vials).

Product Capacity

The purification capacity of the Anti-V5-tag Gel varies depending upon the V5-tagged protein. For example, 1 mL of Anti-V5-tag Gel bound 1.64 mg of a V5-tagged protein (32 kDa) and eluted 1.26 mg of purified protein in our hands.

Storage


Store for up to 1 year from date of receipt at 2-8 °C. Do not freeze.

Material Preparation

Prepare the following reagents before affinity purification.

- Lysis buffer : Suitable Lysis buffer varies with a kind of the V5-tagged protein
(see **Additional Information**).
10-50 mM Tris-HCl (pH 7.5)
100-300 mM NaCl
1-5 mM EDTA
1% NP-40 or Triton X-100
if necessary add Protease Inhibitor Cocktail
(e.g. SIGMA: code P8340, Thermo Fisher Scientific: code 78415)
- Washing buffer : PBS
- Elution buffer : 2 mg/mL V5-tag peptide in PBS.
- Regeneration buffer : 0.17 M Glycine-HCl (pH 2.3)
- Column storage buffer : PBS/Preservative (e.g. 0.1 % ProClin 150 and 0.09% NaN₃)

Procedure Summary

Protocols

- Notes:
1. The Gel is optimized under the native conditions only, and it is not recommended under denaturing conditions and also for purification of aggregated, unstable, and insoluble protein (e.g. inclusion bodies). Proteins solubilized with such as 6 M Guanidine-HCl or 8 M Urea cannot be purified using this Gel (see **Additional Information**).
 2. Cellular debris and particulate matter must be removed prior to purification. The protein extract should be centrifuged (10,000-20,000 × g for 15 minutes) and filtered with a 0.45 µm filter to remove any remaining cells and particulates.
 3. Highly viscous samples containing chromosomal DNA or RNA should be sonicated or treated with nuclease to reduce viscosity.

A. Column preparation

1. Place the empty chromatography column (e.g., Thermo Fisher Scientific, code 29920) on rack or stand.
2. Rinse the column with Washing buffer.
3. Resuspend the Anti-V5-tag Gel by tapping and inverting the vial several times immediately before dispensing. Don't vortex.
4. Transfer the desired volume to the column. Allow the column to drain.
Do not allow the column to dry out.
5. Wash the column with 10 bed volumes of Washing buffer.

B. Column binding

1. Apply the lysate on the top of the column under gravity flow.
Note: The binding efficiency to Anti-V5-tag Gel depends on the conditions such as a kind of the V5-tagged protein, sample flow rate, and temperature. If the coupling efficiency has been low, put a sample through the column several times, or incubate the lysate with the Gel in batch.
2. Collect the flow-through into clean collection tubes.
3. Wash the column with PBS until the OD280 is < 0.01.

C. Elution of the V5-tagged protein

1. Prepare 10 bed volumes of Elution buffer.
2. Remove the bottom cap, and pour off the excess liquid.
3. To exchange the buffer in the column, add 1 bed volume of Elution buffer, and drain the 1 bed volume of the buffer.
4. Cap the bottom, and add 1 bed volumes of Elution buffer so that the Gel does not dry.

5. Allow the column to stand at room temperature for 5 minutes. (As the V5-tag peptide competes with V5-tagged protein on the Gel by this incubation, V5-tagged protein is eluted from the Gel.)
6. Remove the bottom cap, and collect the eluate into clean collection tubes (Fraction 1).
7. Add 1 bed volume of Elution buffer, collect the eluate into clean collection tubes (Fraction 2).
8. Repeat the elution and collection step (Step 6 and 7) another 8 times. (Fraction 1-10 contains V5-tagged Protein.)

Note: Elution efficiency depends on each V5-tagged protein. If the efficiency has been low, please examine numbers of elution, and volume and peptide concentration of Elution buffer.

D. Regeneration and storage

1. Wash the column with 10 bed volumes of Regeneration buffer.
2. Immediately wash the column with 10 bed volumes of Column storage buffer.
3. Store at 2-8 °C in 2 bed volumes of Column storage buffer.


Note: Poured columns containing the anti-V5-tag Gel may be used at least 10 times, depending on the usage conditions.

Related Products:

Please visit our website at <https://ruo.mbl.co.jp/>.

Example of Purification Results

N-terminus V5-tagged β -galactosidase

SDS-PAGE (Coomassie Brilliant Blue Staining)

Human embryonic kidney cells (293T) were transfected with pcDNA-V5-tagged β -galactosidase and cultured for 60 hours. Cells were then lysed in the Lysis buffer (10 mL/100-mm dish \times 5) and purified according to the preceding protocol. Column bed volume was 0.25 mL. Elution was carried out with 1.5 mL of 2 mg/mL V5-tag peptide. Each fraction was 0.25 mL.

Additional Information

Several reagents were examined whether or not they were suitable for use with the V5-tagged Protein PURIFICATION GEL. For example, RIPA buffer could be used for preparation of cell lysate. The results are listed below.

Chaotropic agents

Urea	1 M	Yes
Guanidine-HCl	1 M	Yes

Reducing agents

DTT	10 mM	Yes
2-Mercaptoethanol	10 mM	Yes

Surfactants

Nonionic	Tween-20	5%	Yes
	Triton X-100	5%	Yes
	NP-40	1%	Yes
	Digitonin	1%	Yes
	n-Octyl-β-D-gulcoside	1%	Yes
Zwitterionic	CHAPS	1%	Yes
	CHAPSO	1%	Yes
Anionic	SDS	0.1%	Yes
	Sodium Deoxycholate	0.5%	Yes

Others

NaCl	1 M	Yes
Glycerol	10%	Yes
EDTA	10 mM	Yes

The “Yes” indicates the reagents can be used in the Lysis buffer for this Gel up to the indicated concentration. The “No” indicates the reagents cannot be used in the Lysis buffer for this Gel at the indicated concentration.

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