Revised: September, 2021

For Research Use Only. Not for use in diagnostic procedures.

*For the sensitive detection and quantification of  
mobilized CD107a*

# **IMMUNOCYTO**

## **CD107a Detection Kit**

CODE No. 4844

**MBL** MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.  
A JSR Life Sciences Company  
URL <https://ruo.mbl.co.jp> e-mail [support@mbl.co.jp](mailto:support@mbl.co.jp)

**Before use, thoroughly read these instruction.**

## **Description**

When T-cell and NK-cell recognized target cells, content of the lytic granules are secreted by degranulation. As a result of degradation, CD107 present in the membrane of cytotoxic granules expose onto the cell surface. It follows that the detection of CD107 exposed on the cell surface lead to indirect evaluation of secretion of the lytic granule contents such as perforin, granzymes and granzylsin, which enables rapid assessment of cell-mediated cytotoxicity and detection of antigen specific cytotoxic T-cells.

The method "CD107a mobilization assay" is utilized for its detection. In this method, the cells are cultured in the presence of CD107a monoclonal antibody after antigen stimulation, which enables sensitive detection of CD107a.

IMMUNOCYTO CD107a Detection Kit is a reagent for detection of CD107a exposed on the cell surface after antigen stimulation by CD107a mobilization assay.

## **Intended Use**

IMMUNOCYTO CD107a Detection Kit is intended for the detection of CD107a in research samples.

**For research use only. Not for use in *in vitro* diagnostic procedure for clinical diagnosis.**

## **Product Components**

Each kit contains:

The following reagents are sufficient for 50 assays.

Materials	Quantity
FITC labeled anti-human CD107a monoclonal antibody	240 µL x 1 vial
FITC labeled mouse IgG1 isotype control	240 µL x 1 vial
Monensin	500 µL x 1 vial
Cell suspension solution	25 mL x 2 bottles

## **Precautions**

1. This kit is for research use only. Not recommended or intended for diagnosis of human disease. Not for use in diagnostic procedure.
2. FITC labeled anti-human CD107a monoclonal antibody and FITC labeled isotype control mouse IgG1 contain 0.09% sodium azide as a preservative. Sodium azide has been reported to form potentially explosive lead or copper azide in plumbing. To prevent forming of these compounds, flush drains thoroughly after disposing of solutions containing sodium azide.
3. Monensin contains 70% EtOH. EtOH is flammable.
4. Cell suspension solution contains formaldehyde as fixative. Formaldehyde is harmful, toxic, allergenic and a suspected carcinogen. Avoid inhalation, contact with eyes, skin and clothing.

## Expiration

Please see the label of this kit.

## Storage Conditions

All kit components must be stored at 2-8°C.

## Materials required but not supplied

- Phosphate Buffered Saline (PBS)
- 5 mL test tubes
- Centrifuge
- Aspirator
- Pipettes and disposable pipette tips
- 96-well round microplate
- Medium (RPMI 1640 supplemented with 5-10% serum and 50 IU/mL IL-2)
- T-cell Stimulants (T-cell receptor epitope peptides, OKT3, etc.)
- Flow cytometer

## Protocol

1. Cells are resuspended to final concentration of  $2 \times 10^6$  cells/mL in the medium (RPMI 1640 supplemented with 5-10% serum and 50 IU/mL IL-2).
2. Add the 150  $\mu$ L of cell suspensions to each 4 wells of 96-well round microplate.
3. Add Stimulant, FITC labeled CD107a monoclonal antibody and FITC labeled mouse IgG1 isotype control to each well as the table described below. Monensin prevent reduction of fluorescence of FITC labeled monoclonal antibody which is internalized into an acidic endosomal or lysosomal compartment by neutralization of the pH within endosomes and lysosomes. However CD107a can be detected sufficiently without monensin.

	Cell suspension	Monensin	Stimulant	Medium	CD107a-FITC	IgG1-FITC
well 1	150 $\mu$ L	(2 $\mu$ L)	50 $\mu$ L		2 $\mu$ L	
well 2	150 $\mu$ L	(2 $\mu$ L)	50 $\mu$ L			2 $\mu$ L
well 3	150 $\mu$ L	(2 $\mu$ L)		50 $\mu$ L	2 $\mu$ L	
well 4	150 $\mu$ L	(2 $\mu$ L)		50 $\mu$ L		2 $\mu$ L

\* Stimulant; T-cell receptor epitope peptide (0.04-4  $\mu$ g/mL)

OKT3 (4  $\mu$ g/mL)


etc.

Stimulant should be diluted with Medium.

4. Mix gently by pipetting and Incubate for 1-5 hours in a 5% CO<sub>2</sub> incubator at 37°C.

5. Transfer all cells to test tubes and add 2 mL of PBS (2% FCS and 0.1  $\text{NaN}_3$ ) to each tube.
6. Centrifuge 400 X  $g$  for 5 minutes and Decant supernatant
7. Add 100  $\mu\text{L}$  of PBS (2% FCS and 0.1%  $\text{NaN}_3$ ) to each tube after tapping.
8. If you want to multi-stain the cells with cell surface marker, add the surface marker in this step. For example, add 10  $\mu\text{L}$  of PC5 labeled CD8 and 10  $\mu\text{L}$  of PE labeled MHC tetramer.
9. Mix gently and Incubate all tubes for 15 minutes at room temperature in the dark.
10. Add 2 mL of PBS (2% FCS and 0.1%  $\text{NaN}_3$ ) to each tube.
11. Centrifuge 400 X  $g$  for 5 minutes and decant supernatant.
12. Add 200-500  $\mu\text{L}$  of Cell suspension solution to each tube after tapping, and mix gently.
13. Analysis by flow cytometry.

## Assay Example

**Figure Detection of CD107a exposed on the cell surface after antigen stimulation**  
HLA-A\*2402 restricted CTL line was stimulated by epitope peptide QYDPVAAALF (0.01  $\mu$ g/mL) in the presence of FITC-labeled CD107a monoclonal antibody and cultured for 5 hours at 37°C. After culture, the cells were stained with CD8-PC5 and HLA-A\*2402 CMV pp65 tetramer-PE

Upper column and second column; Not stimulated

Third column and bottom column; Stimulated with pp65 epitope peptide

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