CODE No. 7625

Revised: March, 1 2022 ver. 9

For Research Use Only. Not for use in diagnostic procedures.

Quantitative test kit for Mouse IL-18

Mouse IL-18 ELISA Kit

CODE No. 7625 96 wells

CONTENTS

1. Intended Use -----	1
2. Summary and Explanation-----	1
3. Principle -----	2
4. Materials provided -----	2
5. Storage and Stability -----	2
6. Materials and equipment required -----	3
7. Precautions -----	3
8. Procedure -----	4
9. Performance Characteristics -----	7
10. References -----	9

CODE No. 7625

Before use, thoroughly read these Instructions.

Intended Use

The **Mouse IL-18 ELISA Kit** is based on sandwich ELISA and is capable of measuring mouse IL-18.

For Research Use Only. Not for use in diagnostic procedures.

Summary and Explanation

Interleukin 18 (IL-18) is an 18 kDa cytokine which is identified as a costimulatory factor for production of interferon- γ (IFN- γ) in response to toxic shock. It shares functional similarities with IL-12. IL-18 is synthesized as a precursor 24 kDa molecule without a signal peptide and must be cleaved to produce an active molecule. IL-1 converting enzyme (ICE, caspase-1) cleaves pro-IL-18 at aspartic acid in the P1 position, producing the mature, bioactive peptide that is readily released from the cells. It has been reported that IL-18 is produced from Kupffer cells, activated macrophages, keratinocytes, intestinal epithelial cells, osteoblasts, adrenal cortex cells and murine diencephalons. IFN- γ is produced by activated T and NK cells and plays critical roles in the defense against microbial pathogens. IFN- γ activates macrophages, enhances NK activity and B cell maturation, proliferation and Ig secretion, induces MHC class I and II antigens expression, and inhibits osteoclast activation.

IL-18 acts on T helper 1-type (Th1) cells, and in combination with IL-12 strongly induces production of IFN- γ by these cells. Pleiotropic effects of IL-18 have also been reported, including enhancement production of IFN- γ and GM-CSF in peripheral blood mononuclear cells, production of T helper type 1 cytokines, IL-2, GM-CSF and IFN- γ in T cells, enhancement of Fas ligand expression by Th1 cells.

“Mouse IL-18 ELISA Kit” is the reagent for measuring mouse IL-18 specifically with high sensitivity by ELISA.

1 ng/ml of various cytokines, such as mouse IFN- α , IFN- γ , IL-1 β , IL-4, IL-5, IL-6, IL-10, IL-12, GM-CSF and human IL-18 were not measured by this ELISA. The results were all below the detection limit of 25.0 pg/mL.

Principle

The **Mouse IL-18 ELISA Kit** measures Mouse IL-18 by sandwich ELISA. This assay uses two monoclonal antibodies against two different epitopes of Mouse IL-18. Samples to be measured or standards are incubated in the microwells coated with anti-Mouse IL-18 monoclonal antibody. After washing, the peroxidase conjugated anti-Mouse IL-18 monoclonal antibody is added into the microwells and incubated. After another washing, the substrate reagent is mixed with the chromogen and allowed to incubate for an additional period of time. An acid solution is added to each microwell to terminate the enzyme reaction and to stabilize the color development. The optical density (O.D.) of each microwell is measured at 450 nm using a microplate reader.

The concentration of Mouse IL-18 is calibrated from a dose response curve based on the reference standards.

Materials provided

Each kit contains:

Name	Materials	Quantity (for 1 kit)
Microwells	microwell strips coated with anti-Mouse IL-18 antibody	8-well strip × 12 strips
Mouse IL-18 calibrator	Mouse IL-18 (lyophilized powder)	2 vials
Conjugate reagent	Peroxidase conjugated anti-Mouse IL-18 antibody (× 101)	0.2 mL × 1 vial
Conjugate diluent	Buffer for diluting Conjugate reagent (ready-for-use)	24 mL × 1 bottle
Assay diluent	Buffer for diluting samples* (ready-for-use)	30 mL × 1 bottle
Wash concentrate	Buffer for washing microwells (× 10)	100 mL × 1 bottle
Substrate reagent	TMB/H ₂ O ₂ solution (ready-for-use)	20 mL × 1 bottle
Stop solution	0.5 M H ₂ SO ₄ (ready-for-use)	20 mL × 1 bottle

* Use the Assay diluent as 0 pg/mL of IL-18 standard.

Storage and Stability

All kit components must be stored at 2-8°C. All reagents are stable for 12 months after manufacturing when stored at the indicated conditions.

Please see the label on this kit box.

Materials and equipment required

- Microplate reader (450 nm)
- Automatic plate washer or wash bottle
- Adjustable micropipette
- Multichannel micropipette
- 96-well polyvinyl plate
- Microplate holder
- Uncoated microwell strips (for using an automatic plate washer)
- Disposable reagent vessels
- Paper towels
- Distilled water

Precautions

1. Allow all the components to come to room temperature (20-25°C) before use.
2. **Microwells** that are not immediately required should be returned to the ziplock pouch, which must be carefully resealed to avoid moisture absorption.
3. Fresh samples should be used. Aliquot each sample and store below -20°C if necessary. Avoid repeated freezing and thawing. Never store the samples at 4°C, as samples might be affected by storage of this temperature.
4. When **Wash concentrate** is stored at 2-8°C, some precipitation or turbidity may appear. However it does not affect the reagent efficiency. Allow the **Wash concentrate** to come to room temperature and mix to re-dissolve the precipitate.
5. **Assay diluent** contains sodium azide (0.09%) as preservative. Azide may react with copper or lead in plumbing system to form explosive metal azide. Therefore, always flush with plenty of water into a drain when disposing materials containing azide.
6. **Stop solution** is 0.5 M sulfuric acid. As it is a corrosive material, protect eyes and skin and handle with care.
7. This kit is intended for research use only. Not for use in diagnostic procedures.

Procedure

◆ Preparation of Reagents

1. Wash solution

Prepare “**Wash solution**” by diluting **Wash concentrate, 1:10** with distilled water prior to use. (e.g., add 100 mL of **Wash Concentrate** to 900 mL of distilled water.) This “**Wash solution**” is stable for 2 weeks at 4°C.

2. Conjugate solution

Prepare “**Conjugate solution**” by diluting the **Conjugate reagent, 1:101** with the **Conjugate diluent** prior to use. (e.g., add 10 µL of **Conjugate reagent** to 1,000 µL of **Conjugate diluent**.)

* Prepare only a sufficient amount of the “**Conjugate solution**” for the assay because the diluted **Conjugate reagent** is not stable.

* Use disposable new pipette and vessel to avoid contamination of microbe.

3. Standards

Reconstitute “**Standards**” by dissolving lyophilized **Mouse IL-18 calibrator** with the **Assay diluent**. The volume of **Assay diluent**, and precise procedures are indicated in the attached lot specific document “*Preparation of Standards*”.

* If storage is needed after reconstitution, prepare appropriate aliquots and freeze them below -20°C. Avoid repeated freezing and thawing.

4. Other reagents are ready for use.

◆ Preparation of samples

1. Dilution

1) Dilute each sample with **Assay Diluent**.

e.g. Mouse serum

1:5 with **Assay Diluent**, e.g. by adding 50 µL of sample to 200 µL of **Assay Diluent**. Sample dilution may vary between different specimens. Appropriate sample dilution should be established by each investigator.

Since Mouse IL-18 calibrator contains FCS, any sample which does not contain serum component (such as FCS) may result different reactivity from the sample that contains serum component (e.g. mouse serum, mouse plasma and mouse derived cell culture supernatant containing FCS). Thus, this kit IS NOT SUITABLE FOR MEASURING THE SAMPLE WHICH DOES NOT CONTAIN ANY SERUM COMPONENT.

2) Add 150 µL of prepared samples and **standards** to 96-well polyvinyl plate as the same order of assay run.

2. Storage

Fresh samples should be used. Aliquot each sample into new plastic tube and store below -20°C if necessary. Avoid repeated freezing and thawing.

◆ Assay procedure

Duplicate assay is recommended. A standard curve must be run with each assay.

STEP 1. (Sample incubation)

1) Transfer 100 µL of each sample to Mouse IL-18 antibody coated **Microwells** simultaneously using multichannel pipette.

*Reaction starts on pipetting to Mouse IL-18 antibody coated microwell. Pipetting should be completed as quickly as possible.

2) Cover the plate and incubate for 60 minutes at room temperature (20-25°C).

STEP 2. (Washing)

Aspirate or discard the well contents. Fill the wells with **Wash solution** and then completely aspirate or discard the contents. Wash the wells 4 times with **Wash solution** using washing bottle. When autowasher is used, wash 4 times.

*Each laboratory is recommended to confirm its own appropriate washing times and set-up.

*Washing buffer should be used at room temperature (20-25°C).

STEP 3. (Conjugate incubation)

1) Pour **Conjugate solution** into the vessel. After removing Wash solution remained completely, pipette 100 μL of **Conjugate solution** to each well with multichannel pipette.

*To avoid dry up microwells, the Conjugate solution must be dispensed into the microwells as soon as remove the Wash solution

2) Cover the plate and incubate for 60 minutes at room temperature (20-25°C).

STEP 4. (Washing)

Wash the wells again following the **STEP 2.** procedure.

STEP 5. (Substrate incubation)

1) Pour **Substrate reagent** into the vessel. Add 100 μL of **Substrate reagent** (ready-for-use) to each well.

***Substrate reagent** should be used at room temperature (20-25°C).

*This vessel should be different from the one which was used for pouring **Conjugate solution**.

*Use disposable new pipette and vessel, as **Substrate reagent** is easy to be oxidized by metal ions and be contaminated by microbes.

*If **Substrate reagent** is poured into the vessel from the bottle, do not return to the bottle.

*To avoid dry up microwells, the **Substrate reagent** must be dispensed into the microwells as soon as remove the Wash solution.

2) Cover the plate and incubate for 30 minutes at room temperature (20-25°C).

STEP 6. (Stopping reaction)

Pour **Stop solution** (ready-for-use) into the vessel. Add 100 μL of **Stop solution** to each well with a multichannel pipette.

◆ **Reading**

Read the absorbance of each well at 450 nm using a plate reader. If a dual wavelength plate reader is available, set the test wavelength at 450 nm and the reference at 620 nm.


* Reading should be done within 30 minutes after stopping the reaction.

* Ensure that the back of the plate is clean and dry, and that no air bubbles are present on the surface of the liquid in the wells before reading.

◆ **Calculation of results**

- 1) Calculate the mean absorbance value of each **IL-18 standard**.
 - 2) Plot on the semi-log graph paper and construct a calibration curve.
[Absorbance on the vertical axis, concentration (pg/mL) on the horizontal axis]
 - 3) Read the pg/mL corresponding to the absorbance of sample from the standard curve.
Report the IL-18 concentration of samples by multiplying dilution factor. (e.g., Mouse sera; x 5)
- * If absorbance of sample exceeds the one of the 1,000 pg/mL standard, dilute the sample and measure again.

◆ **Example of standard curve**

Performance Characteristics

◆ **Sensitivity**

The sensitivity of the assay is 25.0 pg/mL.

The minimum detection limit estimated by serial dilution was 25.0 pg/mL since the mean + 2SD of the 12.5 pg/mL was lower than the mean - 2SD of the 25.0 pg/mL.

◆ **Reproducibility**

1. **Intra-assay**

Intra-assay reproducibility was determined by assaying the sera 8 times.

IL-18 concentrations of the serum samples were calculated as described in calculation of results in assay procedure.

Sample	serum 1	serum 2	serum 3	serum 4	serum 5	serum 6
Number of replicates	8	8	8	8	8	8
Mean (pg/ml)	3554.4	3285.9	2791.7	2212.8	992.3	625.3
C.V. (%)	2.5	2.0	5.8	5.3	3.1	3.8

2. **Inter-assay**

Inter-assay reproducibility was determined by 5 independent assays.

IL-18 concentrations of the serum samples were calculated as described in “◆ **Calculation of results**”.


Sample	serum 1	serum 2	serum 3	serum 4	serum 5
Number of determinations	5	5	5	5	5
Mean (pg/ml)*	3500.2	3241.7	2557.9	2066.8	842.6
C.V. (%)	7.4	4.4	3.8	5.4	8.5

* From six replicates of each serum sample in five separate assays.

◆ **Dilution test**

Serum samples were diluted with **Assay diluent**.

IL-18 concentrations of the serum samples were calculated as described in calculation of results in assay procedure.

◆ Recovery test

Recombinant mouse IL-18 was added to samples at different concentrations.

IL-18 concentrations of the serum samples were calculated as described in calculation of results in assay procedure.

serum 1

(A) Additional rmIL-18* (pg/mL)	IL-18 concentration observed (pg/mL)	(B) recovery (pg/mL)	(B/A) recovery (%)
0.0	956.0	-	-
536.0	1497.3	541.3	101.0
960.5	1873.3	917.3	95.5

serum 2

(A) Additional rmIL-18* (pg/mL)	IL-18 concentration observed (pg/mL)	(B) recovery (pg/mL)	(B/A) recovery (%)
0.0	742.0	-	-
536.0	1279.0	537.0	100.2
960.5	1594.0	852.0	88.7

serum 3

(A) Additional rmIL-18* (pg/mL)	IL-18 concentration observed (pg/mL)	(B) recovery (pg/mL)	(B/A) recovery (%)
0.0	1397.3	-	-
536.0	1921.8	524.5	97.9
960.5	2213.0	815.8	84.9

serum 4

(A) Additional rmIL-18* (pg/mL)	IL-18 concentration observed (pg/mL)	(B) recovery (pg/mL)	(B/A) recovery (%)
0.0	340.8	-	-
536.0	884.5	543.8	101.4
960.5	1278.8	938.0	97.7

* rmIL-18 is abbreviation of recombinant mouse IL-18.

References

- 1) Okamura H., et al. Nature 378, 88-91 (1995)
- 2) Ushio S., et al. J. Immunol. 156, 4274-4279 (1996)
- 3) Micallef M., et al. Eur. J. Immunol. 26, 1647-1651 (1996)
- 4) Tao D., et al. Cell Immunol. 173, 230-235 (1998)
- 5) Taniguchi M., et al. J. Immunol. methods 206, 197-113 (1997)

As this product has been used in many researches, these references are a part of such reports.

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