CODE No. 8445

Revised: February, 2022 ver. 3

For Research Use Only. Not for use in diagnostic procedures.

In situ detection kit for Programmed Cell Death

MEBSTAIN

Apoptosis TUNEL Kit Direct

CODE No. 8445

FCM: 66 tests/stain: 40 tests

CONTENTS

1. Intended Use -----	1
2. Summary and Explanation-----	1
3. Principle -----	1
4. Materials provided -----	2
5. Storage -----	2
6. Expiration -----	2
7. Precautions -----	2
8. Procedure -----	2
9. References -----	8

Before use, thoroughly read these Instructions.

Intended Use

The **MEBSTAIN Apoptosis TUNEL Kit Direct** is intended for the detection of apoptotic cells in research samples. It is suitable for histochemistry, cytospin preparation, and flow cytometry. **For Research Use Only. Not for use in diagnostic procedures.**


Summary and Explanation

Programmed cell death is a selective process of physiological cell deletion. It plays a major role in developmental biology and in maintenance of homeostasis in vertebrates^{8), 9), 12), 13)}. For example, during maturation, T cells recognizing self antigens are destroyed by programmed cell death¹¹⁾. Programmed cell death is morphologically known as apoptosis. Apoptosis is accompanied by condensation of cytoplasm, loss of plasma membrane microvilli, condensation and fragmentation of nuclei, and extensive degradation of chromosomal DNA into oligomers of about 180 bp. It has been considered that fragmentation of nuclear DNA is a biochemical hallmark of apoptosis. The TdT-mediated dUTP-biotin nick end labeling method (TUNEL method) developed by Gavrieli et al.¹⁴⁾ has enabled *in situ* visualization of DNA fragmentation at the single cell level and is considered to be a more sensitive method than conventional morphological techniques¹⁰⁾.

The MEBSTAIN Apoptosis TUNEL Kit Direct is an apoptosis detection kit based on the TUNEL method. The main advantage of this rapid and simple procedure is the use of fluorescein-dUTP to label DNA strand breaks. This allows the direct detection of DNA fragmentation.

Principle

In cells in which apoptosis occurs, the chromatin DNA is cut by endonuclease at linker DNA site between nucleosomes. Then, DNA fragments, which are a number of multimers of nucleosomal units exit in nuclei. In TUNEL method, 3'-OH DNA ends generated by DNA fragmentation is nick end labeled with fluorescein-dUTP, mediated by terminal deoxynucleotidyl transferase (TdT). This method allows specific staining.

Materials provided

Each kit contains;

Name	Materials	Quantity (for 1 kit)
Proteinase K (PK)	Proteinase K	0.4 mL × 1 vial
TdT	TdT	0.1 mL × 1 vial
TdT buffer II	Buffer for preparing TdT solution	2.0 mL × 3 vials
FITC-dUTP	FITC conjugated dUTP (ready-for-use)	0.1 mL × 1 vial
TB buffer (×10)	TB buffer (×10)	50 mL × 1 bottle

Storage

All kit components must be stored at -20°C.

Expiration

Please see the label of this kit.

Precautions

1. Allow all the components to come to room temperature (20-25°C) before use.
2. **TdT** contains 90 mM sodium cacodylate and 90 mM potassium cacodylate. **TdT buffer II** contains 100 mM potassium cacodylate. As cacodylic acid is arsenic compound, do not dispose materials containing TdT and TdT buffer II into a drain.
3. **TdT** contains 2-mercaptoethanol (0.05%). Handle with care.
4. This kit is intended for research use only. Not for use in diagnostic procedures.

Procedure

◆ Preparation of Reagents

1. Proteinase K (PK) solution

Prepare “**PK solution**” by diluting **Proteinase K, 1:25** with PBS. Warm this solution to 37°C prior to use.

Example: Dilute 10 µL of Proteinase K to 240 µL of PBS.

2. TdT solution

Prepare “**TdT solution**” by mixing the **TdT buffer II, FITC-dUTP, and TdT** at **18:1:1** prior to use.

Example: Mix 45 µL of TdT buffer II, 2.5 µL of FITC-dUTP, and 2.5 µL of TdT.

* Prepare “**TdT solution**” on the ice.

3. TB solution

Prepare “**TB solution**” by diluting **TB buffer, 1:10** with distilled water prior to use.

Example: Dilute 50 mL of TB buffer to 450 mL of distilled water.

* After thawing TB buffer, it should be stored at 4°C.

◆ for Immunohistochemistry

Materials required but not provided.

- Xylen
- Ethanol (100%, 90%, 80%)
- Phosphate buffered saline (PBS)
- Distilled water
- Mounting medium (90% glycerin, 10% PBS)
- 4% paraformaldehyde [0.1M NaH₂PO₄, pH7.4]
- Propidium Iodide (0.5 mg/mL) *It is for counterstaining.
- Cover slip
- Moist chamber
- Adequate test tube
- Adjustable micropipette

STEP 1. (Preparation of tissue sections)

- 1) Use glass slides pretreated with silane.
- 2) Fix the samples with 4% paraformaldehyde.

Note

4% paraformaldehyde is recommended for sensitive detection. Do not use ethanol and acetone because they may cause non-specific staining.

3-5 µm thick of sections are suitable for staining.

a) Paraffin embedded tissue sections

Fix tissue samples in 4% paraformaldehyde and embedded in paraffin.

b) Frozen tissue sections

Post-Fix cryosections in 4% paraformaldehyde. After the fixation, wash the fixed cryosection with PBS 2 times for 5 minutes each.

When using frozen tissue sections, be careful not to allow it to dry. Proceed directly **STEP 4.**

(DNA nick end labeling).

STEP 2. (Deparaffinization)

Be careful not to allow tissue sections to dry at each step.

- 1) Deparaffinize the sections with Xylene 3 times for 3 minutes each.
- 2) Wash the slides with 100% Ethanol 3 times for 3 minutes each.
- 3) Wash the slides with 90% Ethanol 3 times for 3 minutes each.
- 4) Wash the slides with 80% Ethanol 3 times for 3 minutes each.
- 5) Wash the slides with distilled water 4 times for 2 minutes each.

STEP 3. (Treatment with proteinase K (PK))

This step is not required for frozen tissue sections.

- 1) Pipette 500 μ L of PBS on the section and incubate it for 30 minutes at 37°C. (Or put the slide into a jar filled with PBS which is warmed to 37°C, and incubate it for 15-30 minutes at 37°C.)
- 2) After removing PBS buffer by tapping off, pipette 250 μ L of **PK solution** on the section and incubate it for 30 minutes at 37°C.
- 3) Wash the slides with distilled water at 4 times for 2 minutes each.

Note

In the meantime, prepare **TdT solution**. Preparation should be in the ice.

STEP 4. (DNA nick end labeling)

- 1) Pipette 50 μ L of **TdT buffer II** on the section and incubate it for 5-10 minutes at room temperature.
- 2) After removing **TdT buffer II** by tapping off, pipette 50 μ L of **TdT solution** and incubate it for 60 minutes at 37°C.
- 3) After removing **TdT solution** by tapping off, immerse the slides in **TB solution** in a coplin jar and incubate it for 15 minutes at room temperature.
- 4) Wash the slides with distilled water 4 times for 2 minutes each.

STEP 5. (Counterstaining)

- 1) Pipette 50 μ L of Propidium Iodide (0.5 μ g/mL) on the section and incubate it for 15-20 minutes at 4°C.
- 2) Wash the slides with PBS for 3-5 minutes.

STEP 6. (Mounting and microscopy)

- 1) Mount with mounting medium (90% glycerin, 10% PBS) under coverslip.
- 2) View by fluorescent microscopy.

◆ for Flow Cytometry

It is recommended that the positive and negative controls are set.

Example: as a negative control; peripheral blood mono nuclear cell or cultured cells, such as Raji cells etc. (more than 95% viability).

as a positive control; those above cells which are treated with DNase I (1 µg/mL) for 60 minutes at 37°C. DNase I treatment should be done after permeabilization.

Materials required but not provided.

- Phosphate buffered saline (PBS)
- PBS containing 0.2% BSA
- 4% paraformaldehyde [0.1M NaH₂PO₄, pH7.4]
- 70% Ethanol
- Distilled water
- Adequate test tube
- Adjustable micropipette

STEP 1. (Fixation of cells)

- 1) Wash the cells with PBS containing 0.2% BSA several times.
- 2) Fix the cells with 4% paraformaldehyde for 30 minutes at 4°C.
- 3) Wash the cells with PBS containing 0.2% BSA 2 times.

Note

4% paraformaldehyde is recommended for sensitive detection. Do not use ethanol and acetone because they may cause non-specific staining.

STEP 2. (Permeabilization)

Add 200 µL of 70% Ethanol to the cell pellet, and mix it gently. Then incubate it for 30 minutes at -20°C to permeabilize.

Note

In the meantime, prepare TdT solution. Preparation should be done in the ice.

When the next step can't be done right after this step, cells are suspended in 70% Ethanol and stored at -20°C until test. The test should be done within one week.

STEP 3. (DNA nick end labeling)

- 1) After washing with PBS containing 0.2% BSA 2 times, add 30 μL of **TdT solution** to the cell pellet and incubate it for 60 minutes at 37°C. As negative control, add 30 μL of **negative TdT solution (TdT buffer II and FITC-dUTP at 19:1 (28.5 μL : 1.5 μL))** to the cell pellet and incubate it for 60 minutes at 37°C.
- 2) Wash the cells with PBS containing 0.2% BSA 2 times.
- 3) Suspend the cell pellet to 200-500 μL of PBS containing 0.2% BSA and analyze by a flow cytometer.

◆ for Cytospin samples

It is recommended that the positive and negative controls are set.

Example: as a negative control; peripheral blood mono nuclear cell or cultured cells, such as Raji cells etc. (more than 95% viability).

as a positive control; those above cells which are treated with DNase I (1 $\mu\text{g}/\text{mL}$) for 60 minutes at 37°C. DNase I treatment should be done after permeabilization.

Materials required but not provided.

- PBS containing 2% fetal calf serum (FCS) and 0.09% NaN_3
- Distilled water containing 0.2% BSA
- PBS-T (0.5% Tween-20 and 0.2% BSA in PBS)
- 4% paraformaldehyde [0.1M NaH_2PO_4 , pH7.4]
- Propidium Iodide (0.5 mg/mL) *It is for counterstaining.
- Mounting medium (90% glycerin, 10% PBS)
- Cover slip
- Moist chamber
- Adequate test tube
- Adjustable micropipette

STEP 1. (Cytospin)

Wash the cells with PBS (2% FCS, 0.09% NaN₃) several times and cytopsin according to the normal method.

Condition of cytopsin is depend on the instrument used.

Number of cells: 50-100 μ L of cell suspension (1×10^5 cells/mL) per slide

Rotation: 300-500 \times g

Time: 5-10 minutes

STEP 2. (Fixation of cells)

- 1) Dry the cytopsin slide for about 60 minutes by airing gently.
- 2) Fix the cells with 4% paraformaldehyde for 15 minutes at 4°C.

Note

4% paraformaldehyde is recommended for sensitive detection. Do not use ethanol and acetone because they may cause non-specific staining.

STEP 3. (Permeabilization)

After removing 4% paraformaldehyde by tapping off, pipette PBS-T (0.5% Tween-20 and 0.2% BSA in PBS) onto the cells, and incubate it for about 15 minutes at room temperature.

Note

In the meantime, prepare TdT solution and negative control. Preparation should be done on the ice.

STEP 4. (DNA nick end labeling)

- 1) Wash the slide with distilled water 3 times.
- 2) After removing distilled water by tapping off, pipette 50 μ L of **TdT solution** onto the cells and incubate it for 60 minutes at 37°C. As negative control, add 50 μ L of **negative TdT solution (TdT buffer II and FITC-dUTP at 19:1 (28.5 μ L: 1.5 μ L))** onto the cells and incubate it for 60 minutes at 37°C.
- 3) After removing TdT solution by tapping off, wash the slides with PBS at 3 times.

STEP 5. (Counterstaining)

- 1) Pipette 50 μ L of Propidium Iodide (0.5 μ g/mL) on the section and incubate it for 15-20 minutes at 4°C.
- 2) Wash the slides with PBS 3 minutes.

STEP 6. (Mounting and microscopy)

- 1) Mount with mounting medium (90% glycerin, 10% PBS) under coverslip.
- 2) View by fluorescent microscopy.

References

- 1) Kozaki, K., *et al.*, *Cancer Res.* **68**, 2094-2105 (2008)
- 2) Hara, H., *et al.*, *J. Immunol.* **180**, 7859-7868 (2008)
- 3) Numata, K., *et al.*, *J. Immunol.* **178**, 3777-3785 (2007)
- 4) Yano, T., *et al.*, *Mol. Cell. Biol.* **26**, 4474-4488 (2006)
- 5) Pati, D., *et al.*, *Mol. Cell. Biol.* **22**, 8267-8277 (2002)
- 6) Yoshikawa, R., *et al.*, *Cancer Res.* **61**, 1029-1037 (2001)
- 7) Negoro, S., *et al.*, *Circulation* **103**, 555-561 (2001)
- 8) Steller, H. *Science* **267**, 1445-1449 (1995)
- 9) Mori, C., *et al.*, *Anat. Rec.* **242**, 103-110 (1995)
- 10) Carbonari, M., *et al.*, *Cytometry* **22**, 161-167 (1995)
- 11) Krammer, P. H., *et al.*, *Curr. Opin. Immunol.* **6**, 279-289 (1994)
- 12) Mori, C., *et al.*, *Anat. & Embryol.* **190**, 21-28 (1994)
- 13) Wyllie, A. H., *et al.*, *Br. J. Cancer* **67**, 205-208 (1993)
- 14) Gavrieli, Y., *et al.*, *J. Cell Biol.* **119**, 493-501 (1992)

For more information, please visit our web site.

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