



Met Positive Control

(Human, recombinant protein expressed in Sf9)
Cat# CY-E1080

Lot No.
For 100 Assays
100 units (1 unit/ μ L)

Product Description:

Catalytic domain of human Met (HGF receptor), corresponding to 1076-1370 a.a. containing an N-terminal GST tag and a C-terminal His tag, expressed in recombinant baculovirus infected sf9 cells. Purified by sequentially using GSH agarose and Ni-NTA agarose chromatography. The Met Positive control is designed to use for CycLex Met Kinase Assay/Inhibitor Screening Kit (Cat# CY-1080). The Met Positive Control should be added to the well at 1 unit/well. For instance, diluted positive control 1:10, use 10 μ L for 1 assay. Unused Met Positive Control should be stored at -70°C.

Product Size: 100 units/100 μ L

Formulation: The Met Positive Control is supplied frozen in a buffer containing 20 mM Hepes-KOH (pH 7.5), 1 % BSA, 1 mM EDTA, 2 mM DTT, 50 mM NaCl, 0.03 % Brij35 and 50 % glycerol.

Source: Human Met containing N-terminal GST-tag and C-terminal His tag, expressed in sf9 cells.

Molecular Weight: Met Positive Control demonstrates a single 59 kDa bands by SDS-PAGE analysis.

Purity: Met Positive Control is greater than 75 % pure as determined by SDS-PAGE analysis.

Substrates: Met phosphorylates poly[Glu, Tyr] 4:1 as an exogenous substrate.

Inhibitors: PHA-665752 and SU11274 are known as selective small molecule Met inhibitors ^(10, 11).

Unit Definition: One unit is defined as the amount of kinase required to incorporate 1 nmol of phosphate into the Met (autophosphorylation) under oligomerized-activated condition per 60 minute at 30°C.

Assay Conditions: Assay activity of Met in a 50 μ L reaction containing 20 mM Hepes KOH (pH 7.5), 4 mM $MgCl_2$, 2 mM $MnCl_2$, 1 mM DTT, 50 μ M [γ - ^{32}P] ATP (1 μ Ci), and 4 μ g of CycLex-"Tyrosine kinase-binding module". Start the reaction by adding 10 μ L of the enzyme, diluted 10-fold in a buffer containing 20 mM Hepes KOH (pH 7.5), 1 mM DTT, 0.03 % Brij35. Incubate for 60 minutes at 30°C. Terminate the reaction by adding 600 μ L of cold 10 % TCA solution containing 0.2 % Sodium pyrophosphate and stand on ice for 15 min. Filtrate acid insoluble material through GFC filters (Whatman Inc.), wash 4 times with 1 % TCA and rinse filters with ethanol. Dry filters and count in a liquid scintillation counter.

Storage and Stability: Stable for 12 months at -70°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot enzyme to avoid repeated freezing and thawing.



Met Positive Control

Product Data Sheet

For Research Use Only, Not for use in diagnostic procedures



References:

1. Cooper, C.S., et al. Molecular cloning of a new transforming gene from a chemically transformed human cell line. *Nature*. **311**:29–33, 1984
2. Bottaro, D.P., et al. Identification of the hepatocyte growth factor receptor as the c-met proto-oncogene product. *Science*. **251**:802–804, 1991.
3. Schmidt, L., et al. Germline and somatic mutations in the tyrosine kinase domain of the *MET* proto-oncogene in papillary renal carcinomas. *Nat. Genet.* **16**:68–73, 1997.
4. Schmidt, L., et al. Novel mutations of the *MET* proto-oncogene in papillary renal carcinomas. *Oncogene*. **18**:2343–2350, 1999.
5. Olivero, M., et al. Novel mutation in the ATP-binding site of the *MET* oncogene tyrosine kinase in a HPRCC family. *Int. J. Cancer*. **82**:640–643, 1999.
6. Park, W.S., et al. Somatic mutations in the kinase domain of the Met/hepatocyte growth factor receptor gene in childhood hepatocellular carcinomas. *Cancer Res.* **59**:307–310, 1999.
7. Lee, J.H., et al. A novel germ line juxtamembrane Met mutation in human gastric cancer. *Oncogene*. **19**:4947–4953, 2000.
8. Di Renzo, M.F., et al. Somatic mutations of the *MET* oncogene are selected during metastatic spread of human HNSC carcinomas. *Oncogene*. **19**:1547–1555, 2000.
9. Comoglio, P.M., and Trusolino, L. Invasive growth: from development to metastasis. *J. Clin. Invest.* **109**:857–862, 2002.
10. Christensen JG, Schreck R, Burrows J, Kuruganti P, Chan E, Le P, Chen J, Wang X, Ruslim L, Blake R, Lipson KE, Ramphal J, Do S, Cui JJ, Cherrington JM, Mendel DB. A selective small molecule inhibitor of c-Met kinase inhibits c-Met-dependent phenotypes in vitro and exhibits cytoreductive antitumor activity in vivo. *Cancer Res.* **63**(21):7345-55, 2003
11. Sattler M, Pride YB, Ma P, Gramlich JL, Chu SC, Quinnan LA, Shirazian S, Liang C, Podar K, Christensen JG, Salgia R. A novel small molecule met inhibitor induces apoptosis in cells transformed by the oncogenic TPR-MET tyrosine kinase. *Cancer Res.* **63**(17):5462-9, 2003

For more information, please visit our web site.

<https://ruo.mbl.co.jp/>

MANUFACTURED BY

MBL A JSR Life Sciences Company

MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.

URL: <https://ruo.mbl.co.jp>

E-mail: support@mbl.co.jp

CycLex/CircuLex products are supplied for research use only. CycLex/CircuLex products and components thereof may not be resold, modified for resale, or used to manufacture commercial products without prior written approval from MBL. To inquire about licensing for such commercial use, please contact us via email.