Clear Back

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Clear Back

(Human Fc receptor blocking reagent)

Code No.
MTG-001

Quantity
1 mL (50 tests) x 2 vials

BACKGROUND: This reagent can block non-specific bindings of antibodies to cells caused by human FcR for such experiments as FCM, Fluorescent microscopy, etc. Normal goat IgG in this reagent can also block non-specific bindings caused by any goat secondary antibodies. In an experiment of staining cells with T-Select MHC Tetramer, Clear Back may block non-specific bindings caused by macrophages or endocytosis, resulting in clear staining^{1), 2)}.

FORMULATION: Goat sera contain normal human IgG and 0.09% NaN₃.

*Azide may react with copper or lead in plumbing system to form explosive metal azides. Therefore, always flush plenty of water when disposing materials containing azide into drain.

STRAGE: Clear Back is stable for one year from the date of purchase when stored at 4°C~8°C or at -20°C for longer periods. Avoid repeated freezing and thawing.

APPLICATION:

Blocking buffer: Add 20 µL of Clear Back to 1-5 x 10⁵ cells.

QUANTITY: 1 mL (50 tests) / 1vial x 2

NOTE: Upon storage, a white precipitate might form in this reagent. This precipitate does not affect assay performance. The precipitate may be removed centrifugation and the supernatant can use assay.

REFERENCES:

- 1) Feldman, S. P., *et al.*, *Blood* **61**, 815-818 (1983)
- 2) Daniels, M. A. and Jameson, S. C., *J. Exp. Med.* **191**, 335-345 (2000)
- 3) Fleit, H. B. and Kobasiuk, C. D., *J. Leukoc. Biol.* **49**, 556-565 (1991)
- 4) Kobayashi, M., *et al.*, *J. Ovarian Res.* **7**, 48 (2014)
- 5) Honda, T., *et al.*, *Protein Eng. Des. Sel.* **28**, 53-58 (2015)
- 6) Miyashita, T., *et al.*, *Oncol. Lett.* **14**, 5918-5926 (2017)

This reagent is used in the reference number 4)-6).


The descriptions of the following protocols are examples.
Each user should determine the appropriate condition.

PROTOCOLS:

Flow cytometric analysis for floating cells

- 1) Wash the cells 3 times with washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.09% NaN₃].
- 2) Resuspend the cells with washing buffer (5x10⁶ cells/mL).
- 3) Add 50 µL of the cell suspension into each tube, and centrifuge at 500 x g for 1 minute at room temperature (20~25°C). Remove supernatant by careful aspiration.
- 4) Add 20 µL of Clear Back to the cell pellet after tapping. Mix well and incubate for 5 minutes at room temperature.
- 5) Add 40 µL of the primary antibody at the concentration as suggest in the **APPLICATION** diluted in the washing buffer. Mix well and incubate for 30 minutes at room temperature.
- 6) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove the supernatant by careful aspiration.
- 7) Add 30 µL of the secondary antibody. (The concentration of antibody will depend on the conditions.) Mix well and incubate for 15 minutes at room temperature.
- 8) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove the supernatant by careful aspiration.
- 9) Resuspend the cells with 500 µL of the washing buffer and analyze by a flow cytometer.

Effect of Clear Back on THP-1 cell line


Pink: Blocking with Clear Back
Other colors: Without Clear Back

<Concentration of FITC labeled mouse IgG2a>

| | |
|--------|-----------|
| Pink | 5 µg/mL |
| Blue | 0.1 µg/mL |
| Orange | 0.5 µg/mL |
| Green | 5 µg/mL |

*Clear Back blocks non-specific bindings of antibodies caused by FcR on the THP-1 cell surface.


Comparison of blocking performance between CD16/32 and Clear Back

FITC labeled mouse IgG2a (5 µg/mL) staining THP-1 cells after Clear Back blocking (Green line) and after anti-CD16/32 antibodies blocking (Pink line).

*Clear Back blocks non-specific bindings to cells caused by FcR as much as anti-CD16/32 antibodies do. The anti-CD16/32 antibodies are common-used antibodies for blocking to cells of the bindings caused by FcR.

Clear Back blocks non-specific bindings of anti-CD16 antibody to THP-1 cells.


Green : PE labeled isotype control (mouse IgG1)
 Black : PE labeled anti-CD16

*FcR is expressed on THP-1 cell surface, but CD16 is not. This figure shows that Clear Back blocks non-specific bindings of anti-CD16 antibody to FcR on the surface of THP-1 cells.

Flow cytometric analysis for PBMC (Peripheral Blood Mononuclear Cell)

- 1) Prepare PBMCs in accordance with the usual experimental procedure and resuspend the cells with (5x10⁶ cells/mL).
- 2) Add 20 µL of Clear Back to 50 µL the cell suspension (2.5x10⁵ cells) and incubate 5 minutes at room temperature.
- 3) Incubate the cells with primary antibody. (The concentration of antibody will depend on the conditions.)
- 4) Add 1 mL of washing buffer [PBS containing 2% fetal calf serum (FCS) and 0.09% NaN₃], followed by centrifuge at 500 x g for 1 minute.
- 5) Remove supernatant by careful aspiration.
- 6) Incubate the cells with secondary antibody. (The concentration of antibody will depend on the conditions.)
- 7) Add 1 mL of the washing buffer followed by centrifugation at 500 x g for 1 minute at room temperature. Remove supernatant by careful aspiration.
- 8) Resuspend the cells with 500 µL of the washing buffer and analyze by a flow cytometer.

Effect of Clear Back on PBMC

Pink : Blocking with Clear Back
 Green : Without Clear Back

FITC labeled mouse IgG2a (5 µg/mL) staining PBMC after Clear Back blocking (Pink line) and without Clear Back (Green line).

*Non-specific bindings are not detected on lymphocyte (R2), and Clear Back blocks non-specific bindings of the antibodies to mononuclear cells (R4).