

PEPTIDE

V5-tag peptide

Code No.	Quantity	Form
3315-205	2 mg x 5 vials	Lyophilized

DESCRIPTION:

This synthetic peptide including V5-tag sequence has been designed for the elution of V5-tagged Protein Purification Gel. The peptide sequence is GKPIPNPLLGLDST.

MOLECULAR WEIGHT:

1,421 Da

SUPPLIED:

Lyophilized in PBS

Reconstitute in 1.0 mL of distilled water. After reconstitution, the V5-tag peptide concentration will be 2 mg/mL in PBS.

STORAGE:

Lyophilized V5-tag peptide is stable for one year from the date of shipment when stored at 4°C. After reconstitution, prepare appropriate aliquots and freeze them at -20°C. Avoid repeated freeze/thaw cycles.

RELATED PRODUCTS:

Please visit our website at <https://ruo.mbl.co.jp/>.

191217-3