

PEPTIDE

DDDDK-tag peptide

Code No.	Quantity	Form
3325-205	1 mg x 5 vials	Lyophilized

DESCRIPTION:

This synthetic peptide including DDDDK-tag sequence has been designed for the elution of DDDDK-tagged Protein PURIFICATION KIT. The peptide sequence is DYKDDDDK.

MOLECULAR WEIGHT:

1012 Da

SUPPLIED:

Lyophilized in PBS

Reconstitute in 1.0 mL of distilled water. After reconstitution, the DDDDK-tag peptide concentration will be 1mg/mL in PBS.

STORAGE:

Lyophilized DDDDK-tag peptide is stable for one year from the date of shipment when it is stored at 4°C. After reconstitution, prepare appropriate aliquots and freeze them at -20°C. Avoid repeated freeze/thaw cycles.

For more information, please visit our web site. <https://ruo.mbl.co.jp/>

200908-4