

MONOCLONAL ANTIBODY

Anti-Ring1B mAb

Code No.	Clone	Subclass	Quantity	Concentration
D139-3	3-3	Mouse IgG2b	100 μ L	1 mg/mL

BACKGROUND: Polycomb-group (PcG) proteins form multimeric complexes that maintain the state of transcriptional repression of several regulatory genes during development. Ring1B/Rnf2 forms part of a protein complex containing other PcG proteins, such as Mel18, Bmi1, MBLR, MPc3, and the spliceosome protein Sap155, and these complexes associate with chromatin to regulate transcription. Ring1B may also play a role in the regulation of Hox gene expression by PcG complexes. Deletion of Ring1B activity results in gastrulation arrest and cell cycle inhibition.

SOURCE: This antibody was purified from hybridoma (clone 3-3) supernatant using protein A agarose. This hybridoma was established by fusion of mouse myeloma cell P3U1 with Balb/c mouse splenocyte immunized with the recombinant full-length mouse Ring1B.

FORMULATION: 100 μ g IgG in 100 μ L volume of PBS containing 50% glycerol, pH 7.2. No preservative is contained.

STORAGE: This antibody solution is stable for one year from the date of purchase when stored at -20°C.

REACTIVITY: This antibody reacts with human, mouse and hamster Ring1B on Western blotting.

APPLICATIONS:

Western blotting: 1 μ g/mL for chemiluminescence detection system

Immunoprecipitation: 1-5 μ g/200-300 μ L of cell extract

Immunohistochemistry: Not tested*

Immunofluorescence: Not tested*

Immunocytochemistry: Not tested

Flow cytometry: Not tested

Chromatin Immunoprecipitation: Not tested*

*It is reported that this antibody can be used in Immunohistochemistry⁷⁾, Immunofluorescence³⁾ and Chromatin Immunoprecipitation^{1), 2), 4)-6), 8), 9)}.

Detailed procedure is provided in the following **PROTOCOLS**.

SPECIES CROSS REACTIVITY:


Species	Human	Mouse	Hamster
Samples	Jurkat, U937, HPB-ALL	WR19L, NIH/3T3, embryo lysate	CHO, BHK
Reactivity on WB	+	+	+

INTENDED USE:

For Research Use Only. Not for use in diagnostic procedures.

REFERENCES:

- 1) Vallot, C., *et al.*, *Oncogene* **34**, 741-751 (2015) [ChIP]
- 2) Pradeepa, M. M., *et al.*, *Nucleic Acids Res.* **42**, 9021-9032 (2014) [WB, ChIP]
- 3) Casanova, M., *et al.*, *Cell Rep.* **4**, 1156-1167 (2013) [IF]
- 4) Williamson, I., *et al.*, *Development* **139**, 3157-3167 (2012) [ChIP]
- 5) Mendenhall, E. M., *et al.*, *PLoS Genet.* **6**, e1001244 (2010) [ChIP]
- 6) Landeira, D., *et al.*, *Nat. Cell Biol.* **12**, 618-624 (2010) [ChIP]
- 7) Zaaroor-Regev, D., *et al.*, *PNAS* **107**, 6788-6793 (2010) [IP, IHC]
- 8) Barco, R., *et al.*, *PLoS One* **4**, e5060 (2009) [IP, ChIP]
- 9) Wu, X., *et al.*, *Nucleic Acid Res.* **36**, 3590-3599 (2008) [WB, ChIP]
- 10) Sarcinella, E., *et al.*, *Mol. Cell Biol.* **27**, 6457-6468 (2007) [WB]
- 11) Suzuki, M., *et al.*, *Development* **129**, 4171-4183 (2002)
- 12) Atsuta, T., *et al.*, *Hybridoma* **20**, 43-46 (2001)



Western blot analysis of Ring1B expression in several cells using D139-3.

PROTOCOLS:

SDS-PAGE & Western Blotting

- 1) Wash the cells 3 times with PBS and suspend with 10 volume of cold Lysis buffer [50 mM Tris-HCl (pH 7.2), 250 mM NaCl, 0.1% NP-40, 2 mM EDTA, 10% glycerol] containing appropriate protease inhibitors. Incubate it at 4°C with rotating for 30 minutes, then

- sonicate briefly (up to 10 seconds).
- 2) Centrifuge the tube at 12,000 x g for 10 minutes at 4°C and transfer the supernatant to another tube. Measure the protein concentration of the supernatant and add the cold Lysis buffer to make 8 mg/mL solution.
 - 3) Mix the sample with equal volume of Laemmli's sample buffer.
 - 4) Boil the samples for 3 minutes and centrifuge. Load 10 µL of the sample per lane in a 1 mm thick SDS-polyacrylamide gel for electrophoresis.
 - 5) Blot the protein to a polyvinylidene difluoride (PVDF) membrane at 1 mA/cm² for 1 hour in a semi-dry transfer system (Transfer Buffer: 25 mM Tris, 190 mM glycine, 20% MeOH). See the manufacture's manual for precise transfer procedure.
 - 6) To reduce nonspecific binding, soak the membrane in 5% skimmed milk (in PBS, pH 7.2) for 1 hour at room temperature, or overnight at 4°C.
 - 7) Incubate the membrane with primary antibody diluted with PBS, pH 7.2 containing 5% skimmed milk as suggest in the **APPLICATIONS** for 1 hour at room temperature. (The concentration of antibody will depend on condition.)
 - 8) Wash the membrane with PBS-T [0.05% Tween-20 in PBS] (5 minutes x 6 times).
 - 9) Incubate the membrane with the 1:10,000 Anti-IgG (Mouse) pAb-HRP (MBL; code no. 330) diluted with 5% skimmed milk (in PBS, pH 7.2) for 1 hour at room temperature.
 - 10) Wash the membrane with PBS-T (10 minutes x 3 times).
 - 11) Wipe excess buffer on the membrane, then incubate it with appropriate chemiluminescence reagent for 1 minute.
 - 12) Remove extra reagent from the membrane by dabbing with paper towel, and seal it in plastic wrap.
 - 13) Expose to an X-ray film in a dark room for 2 minutes.
 - 14) Develop the film as usual. The condition for exposure and development may vary.


(Positive controls for Western blotting; Mouse embryo, Jurkat, U937, HPB-ALL, WR19L, NIH/3T3, CHO and BHK)

Immunoprecipitation

- 1) Wash the cells 3 times with PBS and suspend with 10 volume of cold Lysis buffer [50 mM Tris-HCl (pH 7.2), 250 mM NaCl, 0.1% NP-40, 2 mM EDTA, 10% glycerol] containing appropriate protease inhibitors. Incubate it at 4°C with rotating for 30 minutes, then sonicate briefly (up to 10 seconds).
- 2) Centrifuge the tube at 12,000 x g for 10 minutes at 4°C and transfer the supernatant to another tube.
- 3) Add 1-5 µg of the anti-Ring1B monoclonal antibody into 250 µL of the supernatant. Mix well and incubate with gentle agitation for 30-120 minutes at 4°C. Add 20 µL of 50% protein A agarose beads resuspended in the Lysis buffer. Mix well and incubate with gentle agitation for 60 minutes at 4°C.

- 4) Wash the beads 3-5 times with ice-cold Lysis buffer (centrifuge the tube at 2,500 x g for 10 seconds).
- 5) Resuspend the beads in 20 µL of Laemmli's sample buffer, boil for 3-5 minutes, and centrifuge for 5 minutes. Use 10 µL/lane for the SDS-PAGE analysis. (See **SDS-PAGE & Western blotting**.)

(Positive control for Immunoprecipitation; U937)



Immunoprecipitation of Ring1B from U937 cells with D139-3 (1) or normal mouse IgG (2). After immunoprecipitated with the antibody, immunocomplex was resolved on SDS-PAGE and immunoblotted with D139-3.

RELATED PRODUCTS:

- | | |
|--------|--|
| D138-3 | Anti-Sap155 mAb (1A5) |
| D221-3 | Anti-Sap155 mAb (16) |
| PD043 | Anti-Phospho-SF3B1 (Sap155) (Ser129) pAb |
| D210-3 | Anti-Histone H2A mAb (C10037) |