

T-Select PEPTIDE

# H-2K<sup>b</sup> MuLV p15E

## Peptide

## KSPWFTTL

Code No.	Quantity	Concentration
TS-M507-P	100 µL	10 mg/mL

For Research Use Only. Not for administration to humans, and not for diagnostic or therapeutic use.

**Peptide Sequence:** KSPWFTTL (8 aa)

**Origin:** MuLV gp70 p15E, 604-611 aa

**MHC Restriction:** H-2K<sup>b</sup>

**Molecular Weight:** 979.14 Da

**Purity:** ≥90%

**Formulation:** 10 mg/mL peptide solution in DMSO  
Handle using aseptic techniques to avoid contamination.

### Storage Conditions

Store at -20°C. Repeated freeze-thaw cycles should be avoided.

### Background

T lymphocytes play a central role in the immune system. Major histocompatibility complex (MHC) is important for T cell recognition of surface antigen. In humans, MHC is also called human leukocyte antigen (HLA). MHCs present epitope peptides derived from both self and non-self protein, and these MHC/peptide complexes are recognized by T cells via their T-cell receptors (TCRs). In general, CD8<sup>+</sup> cytotoxic T cells (CTLs) recognize the MHC class I/peptide complex and CD4<sup>+</sup> helper T cells recognize MHC class II/peptide complex. Epitope peptides binding to class I molecules typically consist of approximately 8 to 10 consecutive amino acid residues. Meanwhile, epitope peptides binding to class II molecules are not constrained in size and can vary from 11 to 30 amino acids long.

T-Select PEPTIDE is a ready-to-use epitope peptide in liquid form and this peptide is used to produce MHC Tetramer reagent. Epitope peptides enable to stimulate CTLs in an antigen-specific manner, and also to induce T-helper immunity. Therefore, T-Select PEPTIDE series are available for T cell stimulation and expansion to examine T cell functions such as cytokine production, expression of cell surface marker, cytotoxic activity, etc.

Please refer to the data sheet of MHC Tetramer reagent which comprises this epitope peptide and H-2K<sup>b</sup> (MBL, PN TS-M507-1) for more details.

### H-2K<sup>b</sup> MuLV p15E Tetramers

TS-M507-1 H-2K<sup>b</sup> MuLV p15E Tetramer-KSPWFTTL-PE  
TS-M507-2 H-2K<sup>b</sup> MuLV p15E Tetramer-KSPWFTTL-APC

### Mouse Strain H-2K Haplotypes

H-2K allele	H-2K <sup>b</sup>	H-2K <sup>d</sup>	H-2K <sup>k</sup>
Mouse strains	C57BL/-, BXSB/Mp, 129/-	BALB/c, DBA/2, NOD	C3H/He, AKR/J

### References for this epitope

- 1) Sijts AJAM, *et al. J Immunol* **152**: 106-116 (1994)
- 2) Smith PM, *et al. Cell Immunol* **175**: 189-198 (1997)
- 3) Rich RF, *et al. J Virol* **73**: 3826-3834 (1999)
- 4) Yang JC and Perry-Lalley D, *J Immunother* **23**: 177-183 (2000)
- 5) Rich RF, *et al. J Immunol* **168**: 2751-2758 (2002)
- 6) Bronte V, *et al. J Immunol* **171**: 6396-6405 (2003)
- 7) Kudo-Saito C, *et al. Clin Cancer Res* **11**: 2416-2426 (2005)
- 8) Kudo-Saito C, *et al. Clin Cancer Res* **11**: 4533-4544 (2005)
- 9) Rich RF, *et al. Virology* **346**: 287-300 (2006)
- 10) Kudo-Saito C, *et al. Cancer Cell* **15**: 195-206 (2009)
- 11) Dubrot J, *et al. Int J Cancer* **128**: 105-118 (2011)
- 12) Vanrell L, *et al. Mol Ther* **19**: 1245-1253 (2011)
- 13) Quetglas JI, *et al. J Immunol* **190**: 2994-3004 (2013)
- 14) Cash H, *et al. Oncotarget* **6**: 36400-36417 (2015)
- 15) Melero I, *et al. Cancer Res* **75**: 497-507 (2015)

### Related Products

Please check our web site (<https://ruo.mbl.co.jp>) for up-to-date information on products and custom MHC Tetramers.

210115-2