



Phospho-Syntide-2 Mouse Monoclonal Antibody

Product Data Sheet

For Research Use Only, Not for use in diagnostic procedures

Phospho-Syntide-2 Mouse Monoclonal Antibody (Clone MS-6E6)

Cat# CY-M1023

100 µg (1 mg/ml x 100 µL)

Clone Name	Applications	Species Cross-reactivity	Molecular Wt.	Source Isotype
MS-6E6	ELISA	N/A	N/A	Mouse IgG2b

Background: Syntide-2, a peptide based on phosphorylation site two of glycogen synthase, is an exogenous substrate designed for CaM-Kinase (calcium/calmodulin dependent kinases) in the serine-threonine kinase family. CaM-kinase II is a multifunctional calcium/calmodulin dependent protein kinase involved in neuronal functions. This anti-Phospho-Syntide-2 monoclonal antibody has been validated with CaM-kinase II, however it has the potential for use in evaluating other serine threonine kinases such as CaM-Kinase IIa, PKC μ , Akt1, Akt2, Akt3, and PKA.

The relative V_{max}/K_m ratios of the known Ca^{2+} -dependent protein kinases for syntide-2 were determined to be as follows: protein kinase II, 100; protein kinase C, 22; phosphorylase kinase, 2; myosin light chain kinase, 0.005.

Specificity/Sensitivity: Phospho-Syntide-2 Antibody detects phosphorylated Syntide-2 only when phosphorylated at serine residue, by ELISA.

Source/Purification: Monoclonal antibody is produced by immunizing mice with a synthetic phosphopeptide Syntide-2, PLARTL(pS)VAGLPGKK, which is synthetic substrate for CaM-Kinase II (Calmodulin dependent protein kinase II). IgG is purified by protein A-sepharose chromatography.

Recommended Antibody Dilutions: ELISA for detection of CaM-Kinase II activity: 1 µg/mL

Storage: Supplied in 20 mM phosphate buffer (pH 7.5), 300 mM NaCl, 50 % glycerol. Store at -20°C.

Applications Key: WB:Western Blotting IP:Immunoprecipitation IHC:Immunohistochemistry IC:Immunocytochemistry F:Flow cytometry E:ELISA FP:Fluorescence Polarization assay


Species Cross-Reactivity Key: H:human M:mouse R:rat Hm:hamster Mk:monkey Mi:mink C:chicken X:*Xenopus* Z:zebra fish All:all species expected Species enclosed in parentheses are predicted to react based on 100% sequence homology. N/A: Not Applicable



Selected Application References:

1. Hashimoto Y, Soderling TR. Calcium . calmodulin-dependent protein kinase II and calcium phospholipid-dependent protein kinase activities in rat tissues assayed with a synthetic peptide. Arch Biochem Biophys. 252(2):418-25, 1987
2. CM Schworer, RJ Colbran, and TR Soderling Reversible generation of a Ca²⁺-independent form of Ca²⁺(calmodulin)- dependent protein kinase II by an autophosphorylation mechanism J. Biol. Chem., 261: 8581 – 8584, 1986

Fig.1 ELISA for measurement of recombinant CaM-kinase II activity using Phospho-Syntide-2 Monoclonal antibody (MS-6E6) in CycLex CaM-kinase II Assay/Inhibitor Screening Kit (Cat# CY-1173)



Related Products

* CaM kinase II Assay/Inhibitor Screening Kit: Cat# CY-1173

* CaM kinase II Positive Control: Cat# CY-E1173

PRODUCED BY

CycLex Co., Ltd.
1063-103 Terasawaoka
Ina, Nagano 396-0002
Japan
Fax: +81-265-76-7618
e-mail: info@cyclex.co.jp
URL: <http://www.cyclex.co.jp>

CycLex/CircuLex products are supplied for research use only. CycLex/CircuLex products and components thereof may not be resold, modified for resale, or used to manufacture commercial products without prior written approval from CycLex Co., Ltd.. To inquire about licensing for such commercial use, please contact us via email.