

Announcement: Price revision

Price Revision of MBL products, 68 items

Dear Customers,

First of all, please accept our sincere gratitude for your continued support.

Due to the recent drastic increase in raw materials costs, we must regretfully inform about raising the prices of MBL products. The new prices go into effect from April 1st, 2020. Allow us to apologize for any inconvenience this may have caused you.

We appreciate your understanding and cooperation.

Sincerely,

MBL

March, 2020

[Effective date]

April 1st, 2020

[Product information]

Please check Appendix

<Appendix>

Code No.	Product	Size	Previous Price	New Price
D278-6	Anti-GP2 (Glycoprotein 2) (Mouse) mAb-Biotin	100 µg/100 µL	¥36,000	¥58,000
K0218-3	Anti-MICA/B (Human) mAb	100 µg/100 µL	¥30,000	¥58,000
M075-11	Mouse IgG1 (isotype control)-Magnetic Beads	20 tests (Slurry: 1 mL)	¥20,000	¥35,000
M075-12	Mouse IgG1 (isotype control)-ALP	50 µL	¥24,000	¥30,000
M075-3	Mouse IgG1 (isotype control)	100 µg/100 µL	¥20,000	¥25,000
M075-3M2	Mouse IgG1 (isotype control) (Functional Grade)	500 µg/100 µL	¥45,000	¥60,000
M075-4	Mouse IgG1 (isotype control)-FITC	50 µg/1 mL	¥12,000	¥30,000
M075-5	Mouse IgG1 (isotype control)-PE	1 mL (50 tests)	¥12,000	¥30,000
M075-6	Mouse IgG1 (isotype control)-Biotin	50 µg/50 µL	¥20,000	¥30,000
M075-8	Mouse IgG1 (isotype control)-Agarose	400 µL	¥20,000	¥35,000
M075-A48	Mouse IgG1 (isotype control)-Alexa Fluor 488	100 µg	¥16,000	¥30,000
M075-A64	Mouse IgG1 (isotype control)-Alexa Fluor 647	100 µg/100 µL	¥16,000	¥30,000
M076-11	Mouse IgG2a (isotype control)-Magnetic Beads	20 tests (Slurry: 1 mL)	¥20,000	¥35,000
M076-12	Mouse IgG2a (isotype control)-ALP	50 µL	¥24,000	¥30,000
M076-3	Mouse IgG2a (isotype control)	100 µg/100 µL	¥20,000	¥25,000
M076-3M2	Mouse IgG2a (isotype control) (Functional Grade)	500 µg/100 µL	¥45,000	¥60,000
M076-4	Mouse IgG2a (isotype control)-FITC	1 mL	¥12,000	¥30,000
M076-5	Mouse IgG2a (isotype control)-PE	1 mL (50 tests)	¥12,000	¥30,000
M076-6	Mouse IgG2a (isotype control)-Biotin	50 µg/50 µL	¥20,000	¥30,000
M076-A48	Mouse IgG2a (isotype control)-Alexa Fluor 488	100 µg/100 µL	¥16,000	¥30,000
M076-A64	Mouse IgG2a (isotype control)-Alexa Fluor 647	100 µL: 1 mg/mL	¥16,000	¥30,000
M077-11	Mouse IgG2b (isotype control)-Magnetic Beads	20 tests (Slurry: 1 mL)	¥20,000	¥35,000
M077-12	Mouse IgG2b (isotype control)-ALP	50 µL	¥24,000	¥30,000
M077-3	Mouse IgG2b (isotype control)	100 µg/100 µL	¥20,000	¥25,000
M077-3M2	Mouse IgG2b (isotype control) (Functional Grade)	500 µg/100 µL	¥45,000	¥60,000
M077-4	Mouse IgG2b (isotype control)-FITC	50 µg/1 mL	¥12,000	¥30,000
M077-5	Mouse IgG2b (isotype control)-PE	1 mL (50 tests)	¥12,000	¥30,000
M077-6	Mouse IgG2b (isotype control)-Biotin	50 µg/50 µL	¥20,000	¥30,000
M077-A48	Mouse IgG2b (isotype control)-Alexa Fluor 488	100 µg/100 µL	¥16,000	¥30,000
M077-A64	Mouse IgG2b (isotype control)-Alexa Fluor 647	100 µg/100 µL	¥16,000	¥30,000
M078-3	Mouse IgG3 (isotype control)	100 µg/100 µL	¥20,000	¥25,000
M078-3M2	Mouse IgG3 (isotype control) (Functional Grade)	500 µg/100 µL	¥45,000	¥60,000

Code No.	Product	Size	Previous Price	New Price
M078-4	Mouse IgG3 (isotype control)-FITC	50 µg	¥12,000	¥30,000
M078-6	Mouse IgG3 (isotype control)-Biotin	50 µg/50 µL	¥20,000	¥30,000
M079-3	Mouse IgM (isotype control)	100 µg/100 µL	¥20,000	¥25,000
M080-3	Rat IgG1 (isotype control)	100 µg/100 µL	¥10,000	¥25,000
M080-3M2	Rat IgG1 (isotype control) (Functional Grade)	500 µg/100 µL	¥45,000	¥60,000
M080-4	Rat IgG1 (isotype control)-FITC	50 µg	¥12,000	¥30,000
M080-5	Rat IgG1 (isotype control)-PE	1 mL(50tests)	¥12,000	¥30,000
M080-A48	Rat IgG1 (isotype control)-Alexa Fluor 488	100 µg	¥16,000	¥30,000
M080-A64	Rat IgG1 (isotype control)-Alexa Fluor 647	100 µg/100 µL	¥16,000	¥30,000
M081-11	Rat IgG2a (isotype control)-Magnetic Beads	20 tests (Slurry: 1 mL)	¥20,000	¥35,000
M081-3	Rat IgG2a (isotype control)	100 µg/100 µL	¥10,000	¥25,000
M081-3M2	Rat IgG2a (isotype control) (Functional Grade)	500 µg/100 µL	¥45,000	¥60,000
M081-4	Rat IgG2a (isotype control)-FITC	50 µg/1 mL	¥12,000	¥30,000
M081-5	Rat IgG2a (isotype control)-PE	1 mL (50 tests)	¥12,000	¥30,000
M081-8	Rat IgG2a (isotype control)-Agarose	20 tests (Gel: 200 µL)	¥20,000	¥35,000
M081-A48	Rat IgG2a (isotype control)-Alexa Fluor 488	100 µg/100 µL	¥16,000	¥30,000
M081-A64	Rat IgG2a (isotype control)-Alexa Fluor 647	100 µg/100 µL	¥16,000	¥30,000
M082-3	Rat IgG2c (isotype control)	100 µg/100 µL	¥10,000	¥25,000
M082-3M2	Rat IgG2c (isotype control) (Functional Grade)	500 µg/100 µL	¥45,000	¥60,000
M082-4	Rat IgG2c (isotype control)-FITC	50 µg	¥12,000	¥30,000
M090-3	Rat IgG2b (isotype control)	100 µg/100 µL	¥10,000	¥25,000
M090-3M2	Rat IgG2b (isotype control) (Functional Grade)	500 µg/100 µL	¥45,000	¥60,000
M090-4	Rat IgG2b (isotype control)-FITC	1 mL	¥12,000	¥30,000
M090-5	Rat IgG2b (isotype control)-PE	1 mL(50tests)	¥16,000	¥30,000
M090-A48	Rat IgG2b (isotype control)-Alexa Fluor 488	100 µg	¥16,000	¥30,000
M090-A64	Rat IgG2b (isotype control)-Alexa Fluor 647	100 µg/100 µL	¥16,000	¥30,000
M189-3	Syrian Hamster IgG (isotype control)	100 µg/100 µL	¥20,000	¥25,000
M199-3	Armenian Hamster IgG (isotype control)	100 µg/100 µL	¥20,000	¥25,000
PM035	Normal Rabbit IgG	500 µg/100 µL	¥10,000	¥25,000
PM035-8	Normal Rabbit IgG-Agarose	Gel: 200 µL	¥15,000	¥35,000
PM067	Normal Guinea Pig IgG	100 µg/100 µL	¥10,000	¥25,000
PM084	Normal Chicken IgY	200 µg/200 µL	¥10,000	¥25,000
PM084-4	Normal Chicken IgY-FITC	100 µg/100 µL	¥12,000	¥30,000
PM084-5	Normal Chicken IgY-PE	1 mL (50 tests)	¥12,000	¥30,000
PM089	Normal Sheep IgG	500 µg/100 µL	¥10,000	¥25,000
PM094	Normal Goat IgG	500 µg/100 µL	¥10,000	¥25,000