

Announcement: Change of vial tube

Change of vial tube for Recombinant Human IL-18 (without BSA) products

Dear Customers,

First of all, we appreciate for your continued support.

Currently, a packaging of our Recombinant Human IL-18 (without BSA) products is Glass vial. With future versions of the vial tube will be plastic vial tube. The material of these vial tube is polypropylene. Also recombinant protein itself is not changed. We will use the new vial tube sequentially from June 2020. The following tables summarize changes of the vial tubes.

Thank you very much for your kind understanding and cooperation.

Sincerely,
MBL
June, 2020

[Packaging changes Description]




Please see Appendix

[Contact]

About Inventory and Price; TEL : 052-238-1960 E-mail : mbl-hanbaikanri@mbl.co.jp

About Technical Support; TEL : 03-6854-3614 E-mail : support@mbl.co.jp

<Appendix>

Products	Code No.: B003-5 Name: Recombinant Human IL-18 (without BSA) Size: 25 µg	Code No.: B003-2 Name: Recombinant Human IL-18 (without BSA) Size: 200 µg
Current Packaging	Vial: Glass 	Vial: Glass 
New Packaging	Lot. 091~ Tube: Polypropylene (PP) Cap: Polypropylene (PP) with O-ring, Ethylene propylene diene rubber (EPDM) 	Lot. 084~ Tube: Polypropylene (PP) Cap: HD-PE (High density polyethylene) 