

Announcement

## Price and code number change, and new release of MHC Tetramer products (502 items)

Dear Customers,

First of all, please accept our sincere gratitude for your continued support.

We would like to notify you that the price and the code number of MHC Tetramer products listed in Appendix 1 have been changed. This change is due to the relocation of the manufacturing site from the United States to Japan or P. R. China. In addition, the composition of some products will change. Please also refer to Appendix 2 for the newly launched products.

Allow us to apologize for any inconvenience caused by the changes and appreciate your understanding.

Sincerely,

MBL

December, 2024

- **Effective date:** December 2, 2024
- **Product changes:** Please see Appendix 1
- **New products:** Please see Appendix 2

### [Inquiries about the products]

(Technical support group) Email: [support@mbi.co.jp](mailto:support@mbi.co.jp)

## Appendix 1: Product changes

### ■ Manufacturer and code number

Changing the manufacturing site from the United States (MBL International) to Japan (MBL) or P. R. China (MBL Beijing Biotech Co., Ltd., a group company), the product code, “TB-“, will be changed to the following.

- Manufactured in Japan: TS- code
- Manufactured in China: TC- code

### ■ Composition of TC- code items

For products with code TC-, the preservative is changed from sodium azide to ProClin 150 as follow.

Current: 10 mM Tris-HCl (pH 8.0), 150 mM NaCl, 0.5 mM EDTA, 0.2% BSA, <0.1% NaN<sub>3</sub>

New: 10 mM Tris-HCl (pH 8.0), 150 mM NaCl, 0.5 mM EDTA, 0.2% BSA, 0.1% ProClin 150

Code No.		Product Name	Price (w/o tax)	
Current	New		Current	New
TB-0101-1	TC-0101-1	HLA-B*15:01 EBV EBNA4 Tetramer-GQGGSPAM-PE	¥190,000	¥209,000
TB-0101-2	TC-0101-2	HLA-B*15:01 EBV EBNA4 Tetramer-GQGGSPAM-APC	¥190,000	¥209,000
TB-0101-4	TC-0101-4	HLA-B*15:01 EBV EBNA4 Tetramer-GQGGSPAM-BV421	¥190,000	¥209,000
TB-0102-1	TC-0102-1	HLA-A*02:01 Insulin B Tetramer-HLVEALYLV-PE	¥190,000	¥209,000
TB-0102-2	TC-0102-2	HLA-A*02:01 Insulin B Tetramer-HLVEALYLV-APC	¥190,000	¥209,000
TB-0102-4	TC-0102-4	HLA-A*02:01 Insulin B Tetramer-HLVEALYLV-BV421	¥190,000	¥209,000
TB-0103-1	TC-0103-1	HLA-A*03:01 Influenza NP Tetramer-ILRGSVAHK-PE	¥190,000	¥209,000
TB-0103-2	TC-0103-2	HLA-A*03:01 Influenza NP Tetramer-ILRGSVAHK-APC	¥190,000	¥209,000
TB-0103-4	TC-0103-4	HLA-A*03:01 Influenza NP Tetramer-ILRGSVAHK-BV421	¥190,000	¥209,000
TB-0105-1	TC-0105-1	HLA-A*03:01 hTERT Tetramer-KLFGVLRK-PE	¥190,000	¥209,000
TB-0105-2	TC-0105-2	HLA-A*03:01 hTERT Tetramer-KLFGVLRK-APC	¥190,000	¥209,000
TB-0105-4	TC-0105-4	HLA-A*03:01 hTERT Tetramer-KLFGVLRK-BV421	¥190,000	¥209,000
TB-0107-1	TC-0107-1	HLA-A*02:01 IGRP Tetramer-LNIDLLWSV-PE	¥190,000	¥209,000
TB-0107-2	TC-0107-2	HLA-A*02:01 IGRP Tetramer-LNIDLLWSV-APC	¥190,000	¥209,000
TB-0107-4	TC-0107-4	HLA-A*02:01 IGRP Tetramer-LNIDLLWSV-BV421	¥190,000	¥209,000
TB-0108-1	TC-0108-1	HLA-A*03:01 Survivin 18-27 (modK) Tetramer-RISTFKNWPK-PE	¥190,000	¥209,000
TB-0108-2	TC-0108-2	HLA-A*03:01 Survivin 18-27 (modK) Tetramer-RISTFKNWPK-APC	¥190,000	¥209,000
TB-0108-4	TC-0108-4	HLA-A*03:01 Survivin 18-27 (modK) Tetramer-RISTFKNWPK-BV421	¥190,000	¥209,000
TB-0109-1	TC-0109-1	HLA-A*03:01 HIV p17 gag Tetramer-RLRPGGKKK-PE	¥190,000	¥209,000
TB-0109-2	TC-0109-2	HLA-A*03:01 HIV p17 gag Tetramer-RLRPGGKKK-APC	¥190,000	¥209,000
TB-0109-4	TC-0109-4	HLA-A*03:01 HIV p17 gag Tetramer-RLRPGGKKK-BV421	¥190,000	¥209,000
TB-0110-1	TC-0110-1	HLA-A*02:01 Mesothelin Tetramer-SLLFLLFSL-PE	¥190,000	¥209,000
TB-0110-2	TC-0110-2	HLA-A*02:01 Mesothelin Tetramer-SLLFLLFSL-APC	¥190,000	¥209,000
TB-0110-4	TC-0110-4	HLA-A*02:01 Mesothelin Tetramer-SLLFLLFSL-BV421	¥190,000	¥209,000
TB-0112-1	TC-0112-1	HLA-A*02:01 Mesothelin Tetramer-VLPLTVAEV-PE	¥190,000	¥209,000
TB-0112-2	TC-0112-2	HLA-A*02:01 Mesothelin Tetramer-VLPLTVAEV-APC	¥190,000	¥209,000
TB-0112-4	TC-0112-4	HLA-A*02:01 Mesothelin Tetramer-VLPLTVAEV-BV421	¥190,000	¥209,000
TB-0113-1	TC-0113-1	HLA-A*02:01 hTERT Tetramer-YLFFYRKS-PE	¥190,000	¥209,000

TB-0113-2	TC-0113-2	HLA-A*02:01 hTERT Tetramer-YLFFYRKSIV-APC	¥190,000	¥209,000
TB-0113-4	TC-0113-4	HLA-A*02:01 hTERT Tetramer-YLFFYRKSIV-BV421	¥190,000	¥209,000
TB-0114-1	TC-0114-1	HLA-A*02:01 hTERT Tetramer-YLQVNSLQTV-PE	¥190,000	¥209,000
TB-0114-2	TC-0114-2	HLA-A*02:01 hTERT Tetramer-YLQVNSLQTV-APC	¥190,000	¥209,000
TB-0114-4	TC-0114-4	HLA-A*02:01 hTERT Tetramer-YLQVNSLQTV-BV421	¥190,000	¥209,000
TB-0116-1	TC-0116-1	HLA-A*02:01 LCMV GPC 447-455 Tetramer-YLVSIFLHL-PE	¥190,000	¥209,000
TB-0116-2	TC-0116-2	HLA-A*02:01 LCMV GPC 447-455 Tetramer-YLVSIFLHL-APC	¥190,000	¥209,000
TB-0116-4	TC-0116-4	HLA-A*02:01 LCMV GPC 447-455 Tetramer-YLVSIFLHL-BV421	¥190,000	¥209,000
TB-0117-4	TC-0117-4	HLA-A*24:02 EBV LMP-2 Tetramer-TYGPVFMCL-BV421	¥190,000	¥209,000
TB-0118-1	TC-0118-1	HLA-A*02:01 HCV NS3 1073-1081 Tetramer-CVNGVCWTV-PE	¥190,000	¥209,000
TB-0118-2	TC-0118-2	HLA-A*02:01 HCV NS3 1073-1081 Tetramer CVNGVCWTV-APC	¥190,000	¥209,000
TB-0118-4	TC-0118-4	HLA-A*02:01 HCV NS3 1073-1081 Tetramer-CVNGVCWTV-BV421	¥190,000	¥209,000
TB-0119-1	TC-0119-1	HLA-B*07:02 Tyrosinase Tetramer-LPWHRFLFLL-PE	¥190,000	¥209,000
TB-0119-2	TC-0119-2	HLA-B*07:02 Tyrosinase Tetramer-LPWHRFLFLL-APC	¥190,000	¥209,000
TB-0119-4	TC-0119-4	HLA-B*07:02 Tyrosinase Tetramer-LPWHRFLFLL-BV421	¥190,000	¥209,000
TB-0120-1	TC-0120-1	HLA-A*02:01 gp100 Tetramer-YLEPGPVTV-PE	¥190,000	¥209,000
TB-0120-2	TC-0120-2	HLA-A*02:01 gp100 Tetramer-YLEPGPVTV-APC	¥190,000	¥209,000
TB-0120-4	TC-0120-4	HLA-A*02:01 gp100 Tetramer-YLEPGPVTV-BV421	¥190,000	¥209,000
TB-0121-1	TC-0121-1	HLA-A*02:01 LCMV ND Tetramer-ALPHIIDEV-PE	¥190,000	¥209,000
TB-0121-2	TC-0121-2	HLA-A*02:01 LCMV ND Tetramer-ALPHIIDEV-APC	¥190,000	¥209,000
TB-0121-4	TC-0121-4	HLA-A*02:01 LCMV ND Tetramer-ALPHIIDEV-BV421	¥190,000	¥209,000
TB-0122-1	TC-0122-1	HLA-A*02:01 HBV S protein Tetramer-FLLTRILTI-PE	¥190,000	¥209,000
TB-0122-2	TC-0122-2	HLA-A*02:01 HBV S protein Tetramer-FLLTRILTI-APC	¥190,000	¥209,000
TB-0122-4	TC-0122-4	HLA-A*02:01 HBV S protein Tetramer-FLLTRILTI-BV421	¥190,000	¥209,000
TB-0123-1	TC-0123-1	HLA-B*07:02 EBV EBNA 6 Tetramer-QPRAPIRPI-PE	¥190,000	¥209,000
TB-0123-2	TC-0123-2	HLA-B*07:02 EBV EBNA 6 Tetramer-QPRAPIRPI-APC	¥190,000	¥209,000
TB-0123-4	TC-0123-4	HLA-B*07:02 EBV EBNA 6 Tetramer-QPRAPIRPI-BV421	¥190,000	¥209,000
TB-0125-1	TC-0125-1	HLA-A*01:01 HCV polyprotein Tetramer-ATDALMTGY-PE	¥190,000	¥209,000
TB-0125-2	TC-0125-2	HLA-A*01:01 HCV polyprotein Tetramer-ATDALMTGY-APC	¥190,000	¥209,000
TB-0125-4	TC-0125-4	HLA-A*01:01 HCV polyprotein Tetramer-ATDALMTGY-BV421	¥190,000	¥209,000
TB-0126-1	TC-0126-1	HLA-A*01:01 Tyrosinase Tetramer-KSDICTDEY-PE	¥190,000	¥209,000
TB-0126-2	TC-0126-2	HLA-A*01:01 Tyrosinase Tetramer-KSDICTDEY-APC	¥190,000	¥209,000
TB-0126-4	TC-0126-4	HLA-A*01:01 Tyrosinase Tetramer-KSDICTDEY-BV421	¥190,000	¥209,000
TB-0127-1	TC-0127-1	HLA-A*02:01 PSCA Tetramer-AILALLPAL-PE	¥190,000	¥209,000
TB-0127-2	TC-0127-2	HLA-A*02:01 PSCA Tetramer-AILALLPAL-APC	¥190,000	¥209,000
TB-0127-4	TC-0127-4	HLA-A*02:01 PSCA Tetramer-AILALLPAL-BV421	¥190,000	¥209,000
TB-0128-1	TC-0128-1	HLA-A*02:01 hTERT Tetramer-ALLTSRLRFI-PE	¥190,000	¥209,000
TB-0128-2	TC-0128-2	HLA-A*02:01 hTERT Tetramer-ALLTSRLRFI-APC	¥190,000	¥209,000
TB-0128-4	TC-0128-4	HLA-A*02:01 hTERT Tetramer-ALLTSRLRFI-BV421	¥190,000	¥209,000
TB-0129-1	TC-0129-1	HLA-B*35:01 NY-ESO-1 Tetramer-MPFATPMEA-PE	¥190,000	¥209,000
TB-0129-2	TC-0129-2	HLA-B*35:01 NY-ESO-1 Tetramer-MPFATPMEA-APC	¥190,000	¥209,000
TB-0129-4	TC-0129-4	HLA-B*35:01 NY-ESO-1 Tetramer-MPFATPMEA-BV421	¥190,000	¥209,000
TB-0130-1	TC-0130-1	HLA-A*24:02 HIV env Tetramer-RYLKQDQLL-PE	¥190,000	¥209,000

TB-0130-2	TC-0130-2	HLA-A*24:02 HIV env Tetramer-RYLKDQQLL-APC	¥190,000	¥209,000
TB-0130-4	TC-0130-4	HLA-A*24:02 HIV env Tetramer-RYLKDQQLL-BV421	¥190,000	¥209,000
TB-0131-1	TC-0131-1	HLA-A*02:01 EBV BALF4 Tetramer-FLDKGTYTL-PE	¥190,000	¥209,000
TB-0131-2	TC-0131-2	HLA-A*02:01 EBV BALF4 Tetramer-FLDKGTYTL-APC	¥190,000	¥209,000
TB-0131-4	TC-0131-4	HLA-A*02:01 EBV BALF4 Tetramer-FLDKGTYTL-BV421	¥190,000	¥209,000
TB-0132-1	TC-0132-1	HLA-A*02:01 topo II Tetramer-FLYDDNQRV-PE	¥190,000	¥209,000
TB-0132-2	TC-0132-2	HLA-A*02:01 topo II Tetramer-FLYDDNQRV-APC	¥190,000	¥209,000
TB-0132-4	TC-0132-4	HLA-A*02:01 topo II Tetramer-FLYDDNQRV-BV421	¥190,000	¥209,000
TB-0133-1	TC-0133-1	HLA-A*02:01 HTLV bZIP Tetramer-AVLDGLLSL-PE	¥190,000	¥209,000
TB-0133-2	TC-0133-2	HLA-A*02:01 HTLV bZIP Tetramer-AVLDGLLSL-APC	¥190,000	¥209,000
TB-0133-4	TC-0133-4	HLA-A*02:01 HTLV bZIP Tetramer-AVLDGLLSL-BV421	¥190,000	¥209,000
TB-0134-1	TC-0134-1	HLA-A*02:01 GPC3 Tetramer-FVGEFFTDV-PE	¥190,000	¥209,000
TB-0134-2	TC-0134-2	HLA-A*02:01 GPC3 Tetramer-FVGEFFTDV-APC	¥190,000	¥209,000
TB-0134-4	TC-0134-4	HLA-A*02:01 GPC3 Tetramer-FVGEFFTDV-BV421	¥190,000	¥209,000
TB-0135-1	TC-0135-1	HLA-A*02:01 HIV Vif Tetramer-GLADQLIHL-PE	¥190,000	¥209,000
TB-0135-2	TC-0135-2	HLA-A*02:01 HIV Vif Tetramer-GLADQLIHL-APC	¥190,000	¥209,000
TB-0135-4	TC-0135-4	HLA-A*02:01 HIV Vif Tetramer-GLADQLIHL-BV421	¥190,000	¥209,000
TB-0136-1	TC-0136-1	HLA-A*02:01 p53 Tetramer-GLAPPQHILRV-PE	¥190,000	¥209,000
TB-0136-2	TC-0136-2	HLA-A*02:01 p53 Tetramer-GLAPPQHILRV-APC	¥190,000	¥209,000
TB-0136-4	TC-0136-4	HLA-A*02:01 p53 Tetramer-GLAPPQHILRV-BV421	¥190,000	¥209,000
TB-0137-1	TC-0137-1	HLA-A*03:01 HPV E6 Tetramer-KLCLRFSLK-PE	¥190,000	¥209,000
TB-0137-2	TC-0137-2	HLA-A*03:01 HPV E6 Tetramer-KLCLRFSLK-APC	¥190,000	¥209,000
TB-0137-4	TC-0137-4	HLA-A*03:01 HPV E6 Tetramer-KLCLRFSLK-BV421	¥190,000	¥209,000
TB-0138-1	TC-0138-1	HLA-A*03:01 HIV nef Tetramer-QVPLRPMTYK-PE	¥190,000	¥209,000
TB-0138-2	TC-0138-2	HLA-A*03:01 HIV nef Tetramer-QVPLRPMTYK-APC	¥190,000	¥209,000
TB-0138-4	TC-0138-4	HLA-A*03:01 HIV nef Tetramer-QVPLRPMTYK-BV421	¥190,000	¥209,000
TB-0139-1	TC-0139-1	HLA-A*24:02 PSA Tetramer-CYASGWGSI-PE	¥190,000	¥209,000
TB-0139-2	TC-0139-2	HLA-A*24:02 PSA Tetramer-CYASGWGSI-APC	¥190,000	¥209,000
TB-0139-4	TC-0139-4	HLA-A*24:02 PSA Tetramer-CYASGWGSI-BV421	¥190,000	¥209,000
TB-0140-1	TC-0140-1	HLA-A*24:02 GPC3 Tetramer-EYILSLEEL-PE	¥190,000	¥209,000
TB-0140-2	TC-0140-2	HLA-A*24:02 GPC3 Tetramer-EYILSLEEL-APC	¥190,000	¥209,000
TB-0140-4	TC-0140-4	HLA-A*24:02 GPC3 Tetramer-EYILSLEEL-BV421	¥190,000	¥209,000
TB-0142-1	TC-0142-1	HLA-B*07:02 HIV gag Tetramer-GPGHKARVL-PE	¥190,000	¥209,000
TB-0142-2	TC-0142-2	HLA-B*07:02 HIV gag Tetramer-GPGHKARVL-APC	¥190,000	¥209,000
TB-0142-4	TC-0142-4	HLA-B*07:02 HIV gag Tetramer-GPGHKARVL-BV421	¥190,000	¥209,000
TB-0143-1	TC-0143-1	HLA-B*07:02 HPV E7 Tetramer-KPTLKEYVL-PE	¥190,000	¥209,000
TB-0143-2	TC-0143-2	HLA-B*07:02 HPV E7 Tetramer-KPTLKEYVL-APC	¥190,000	¥209,000
TB-0143-4	TC-0143-4	HLA-B*07:02 HPV E7 Tetramer-KPTLKEYVL-BV421	¥190,000	¥209,000
TB-0144-1	TC-0144-1	HLA-B*07:02 H-Y Tetramer-SPSVDKARAE-PE	¥190,000	¥209,000
TB-0144-2	TC-0144-2	HLA-B*07:02 H-Y Tetramer-SPSVDKARAE-APC	¥190,000	¥209,000
TB-0144-4	TC-0144-4	HLA-B*07:02 H-Y Tetramer-SPSVDKARAE-BV421	¥190,000	¥209,000
TB-0145-1	TC-0145-1	HLA-B*08:01 HCV NS3 Tetramer-HSKKKCDEL-PE	¥190,000	¥209,000
TB-0145-2	TC-0145-2	HLA-B*08:01 HCV NS3 Tetramer-HSKKKCDEL-APC	¥190,000	¥209,000

TB-0145-4	TC-0145-4	HLA-B*08:01 HCV NS3 Tetramer-HSKKKCDEL-BV421	¥190,000	¥209,000
TB-0146-1	TC-0146-1	HLA-A*02:01 EBV LMP1 Tetramer-YLLEMLWRL-PE	¥190,000	¥209,000
TB-0146-2	TC-0146-2	HLA-A*02:01 EBV LMP1 Tetramer-YLLEMLWRL-APC	¥190,000	¥209,000
TB-0146-4	TC-0146-4	HLA-A*02:01 EBV LMP1 Tetramer-YLLEMLWRL-BV421	¥190,000	¥209,000
TB-0147-1	TC-0147-1	HLA-B*08:01 CMV IE1 Tetramer-QIKVRVDMV-PE	¥190,000	¥209,000
TB-0147-2	TC-0147-2	HLA-B*08:01 CMV IE1 Tetramer-QIKVRVDMV-APC	¥190,000	¥209,000
TB-0147-4	TC-0147-4	HLA-B*08:01 CMV IE1 Tetramer-QIKVRVDMV-BV421	¥190,000	¥209,000
TB-0149-1	TC-0149-1	HLA-B*35:01 HIV RT Tetramer-VPLDEDFRKY-PE	¥190,000	¥209,000
TB-0149-2	TC-0149-2	HLA-B*35:01 HIV RT Tetramer-VPLDEDFRKY-APC	¥190,000	¥209,000
TB-0149-4	TC-0149-4	HLA-B*35:01 HIV RT Tetramer-VPLDEDFRKY-BV421	¥190,000	¥209,000
TB-0150-1	TC-0150-1	HLA-A*02:01 hTERT Tetramer-GLLGASVLGL-PE	¥190,000	¥209,000
TB-0150-2	TC-0150-2	HLA-A*02:01 hTERT Tetramer-GLLGASVLGL-APC	¥190,000	¥209,000
TB-0150-4	TC-0150-4	HLA-A*02:01 hTERT Tetramer-GLLGASVLGL-BV421	¥190,000	¥209,000
TB-0151-1	TC-0151-1	HLA-A*02:01 BA46 Tetramer-GLQHWVPEL-PE	¥190,000	¥209,000
TB-0151-2	TC-0151-2	HLA-A*02:01 BA46 Tetramer-GLQHWVPEL-APC	¥190,000	¥209,000
TB-0151-4	TC-0151-4	HLA-A*02:01 BA46 Tetramer-GLQHWVPEL-BV421	¥190,000	¥209,000
TB-0152-1	TC-0152-1	HLA-A*02:01 p53 Tetramer-KLCPVQLWV-PE	¥190,000	¥209,000
TB-0152-2	TC-0152-2	HLA-A*02:01 p53 Tetramer-KLCPVQLWV-APC	¥190,000	¥209,000
TB-0152-4	TC-0152-4	HLA-A*02:01 p53 Tetramer-KLCPVQLWV-BV421	¥190,000	¥209,000
TB-0153-1	TC-0153-1	HLA-A*02:01 MUC1 Tetramer-LLLTVLTVV-PE	¥190,000	¥209,000
TB-0153-2	TC-0153-2	HLA-A*02:01 MUC1 Tetramer-LLLTVLTVV-APC	¥190,000	¥209,000
TB-0153-4	TC-0153-4	HLA-A*02:01 MUC1 Tetramer-LLLTVLTVV-BV421	¥190,000	¥209,000
TB-0154-1	TC-0154-1	HLA-A*02:01 HHV-8 gB Tetramer-LMWYELSKI-PE	¥190,000	¥209,000
TB-0154-2	TC-0154-2	HLA-A*02:01 HHV-8 gB Tetramer-LMWYELSKI-APC	¥190,000	¥209,000
TB-0154-4	TC-0154-4	HLA-A*02:01 HHV-8 gB Tetramer-LMWYELSKI-BV421	¥190,000	¥209,000
TB-0155-1	TC-0155-1	HLA-A*02:01 Survivin Tetramer-LTLGEFLKL-PE	¥190,000	¥209,000
TB-0155-2	TC-0155-2	HLA-A*02:01 Survivin Tetramer-LTLGEFLKL-APC	¥190,000	¥209,000
TB-0155-4	TC-0155-4	HLA-A*02:01 Survivin Tetramer-LTLGEFLKL-BV421	¥190,000	¥209,000
TB-0156-1	TC-0156-1	HLA-A*02:01 LIVIN Tetramer-QLCPCICRAPV-PE	¥190,000	¥209,000
TB-0156-2	TC-0156-2	HLA-A*02:01 LIVIN Tetramer-QLCPCICRAPV-APC	¥190,000	¥209,000
TB-0156-4	TC-0156-4	HLA-A*02:01 LIVIN Tetramer-QLCPCICRAPV-BV421	¥190,000	¥209,000
TB-0157-1	TC-0157-1	HLA-A*02:01 p53 Tetramer-RMPEAAPPV-PE	¥190,000	¥209,000
TB-0157-2	TC-0157-2	HLA-A*02:01 p53 Tetramer-RMPEAAPPV-APC	¥190,000	¥209,000
TB-0157-4	TC-0157-4	HLA-A*02:01 p53 Tetramer-RMPEAAPPV-BV421	¥190,000	¥209,000
TB-0158-1	TC-0158-1	HLA-A*02:01 CD20 Tetramer-SLFLGILSV-PE	¥190,000	¥209,000
TB-0158-2	TC-0158-2	HLA-A*02:01 CD20 Tetramer-SLFLGILSV-APC	¥190,000	¥209,000
TB-0158-4	TC-0158-4	HLA-A*02:01 CD20 Tetramer-SLFLGILSV-BV421	¥190,000	¥209,000
TB-0159-1	TC-0159-1	HLA-A*02:01 p53 Tetramer-SLPPPGTRV-PE	¥190,000	¥209,000
TB-0159-2	TC-0159-2	HLA-A*02:01 p53 Tetramer-SLPPPGTRV-APC	¥190,000	¥209,000
TB-0159-4	TC-0159-4	HLA-A*02:01 p53 Tetramer-SLPPPGTRV-BV421	¥190,000	¥209,000
TB-0161-1	TC-0161-1	HLA-A*02:01 PSMA Tetramer-VLAGGFLL-PE	¥190,000	¥209,000
TB-0161-2	TC-0161-2	HLA-A*02:01 PSMA Tetramer-VLAGGFLL-APC	¥190,000	¥209,000
TB-0161-4	TC-0161-4	HLA-A*02:01 PSMA Tetramer-VLAGGFLL-BV421	¥190,000	¥209,000

TB-0162-1	TC-0162-1	HLA-A*02:01 IGRP Tetramer-VLFGLGFAI-PE	¥190,000	¥209,000
TB-0162-2	TC-0162-2	HLA-A*02:01 IGRP Tetramer-VLFGLGFAI-APC	¥190,000	¥209,000
TB-0162-4	TC-0162-4	HLA-A*02:01 IGRP Tetramer-VLFGLGFAI-BV421	¥190,000	¥209,000
TB-0163-1	TC-0163-1	HLA-A*02:01 p53 Tetramer-YLGSYGFRL-PE	¥190,000	¥209,000
TB-0163-2	TC-0163-2	HLA-A*02:01 p53 Tetramer-YLGSYGFRL-APC	¥190,000	¥209,000
TB-0163-4	TC-0163-4	HLA-A*02:01 p53 Tetramer-YLGSYGFRL-BV421	¥190,000	¥209,000
TB-0164-1	TC-0164-1	HLA-A*02:01 EZH2 Tetramer-YMCSFLFNL-PE	¥190,000	¥209,000
TB-0164-2	TC-0164-2	HLA-A*02:01 EZH2 Tetramer-YMCSFLFNL-APC	¥190,000	¥209,000
TB-0164-4	TC-0164-4	HLA-A*02:01 EZH2 Tetramer-YMCSFLFNL-BV421	¥190,000	¥209,000
TB-0165-4	TC-0165-4	HLA-A*02:01 Influenza M1 Tetramer-ILGFVFTLTV-BV421	¥190,000	¥209,000
TB-0166-1	TC-0166-1	HLA-A*03:01 gp100 Tetramer-ALLAVGATK-PE	¥190,000	¥209,000
TB-0166-2	TC-0166-2	HLA-A*03:01 gp100 Tetramer-ALLAVGATK-APC	¥190,000	¥209,000
TB-0166-4	TC-0166-4	HLA-A*03:01 gp100 Tetramer-ALLAVGATK-BV421	¥190,000	¥209,000
TB-0167-1	TC-0167-1	HLA-A*24:02 LY6K Tetramer-RYCNLEGPPI-PE	¥190,000	¥209,000
TB-0167-2	TC-0167-2	HLA-A*24:02 LY6K Tetramer-RYCNLEGPPI-APC	¥190,000	¥209,000
TB-0167-4	TC-0167-4	HLA-A*24:02 LY6K Tetramer-RYCNLEGPPI-BV421	¥190,000	¥209,000
TB-0168-1	TC-0168-1	HLA-B*35:01 EBV EBNA1 Tetramer-HPVGEADYFEY-PE	¥190,000	¥209,000
TB-0168-2	TC-0168-2	HLA-B*35:01 EBV EBNA1 Tetramer-HPVGEADYFEY-APC	¥190,000	¥209,000
TB-0168-4	TC-0168-4	HLA-B*35:01 EBV EBNA1 Tetramer-HPVGEADYFEY-BV421	¥190,000	¥209,000
TB-0169-1	TC-0169-1	HLA-A*02:01 BCR-ABL Tetramer-GVRGRVEEI-PE	¥190,000	¥209,000
TB-0169-2	TC-0169-2	HLA-A*02:01 BCR-ABL Tetramer-GVRGRVEEI-APC	¥190,000	¥209,000
TB-0169-4	TC-0169-4	HLA-A*02:01 BCR-ABL Tetramer-GVRGRVEEI-BV421	¥190,000	¥209,000
TB-0170-1	TC-0170-1	HLA-A*02:01 IAPP Tetramer-KLQVFLIVL-PE	¥190,000	¥209,000
TB-0170-2	TC-0170-2	HLA-A*02:01 IAPP Tetramer-KLQVFLIVL-APC	¥190,000	¥209,000
TB-0170-4	TC-0170-4	HLA-A*02:01 IAPP Tetramer-KLQVFLIVL-BV421	¥190,000	¥209,000
TB-0173-1	TC-0173-1	HLA-A*02:01 HPV16 E7 Tetramer-MLDLQPETT-PE	¥190,000	¥209,000
TB-0173-2	TC-0173-2	HLA-A*02:01 HPV16 E7 Tetramer-MLDLQPETT-APC	¥190,000	¥209,000
TB-0173-4	TC-0173-4	HLA-A*02:01 HPV16 E7 Tetramer-MLDLQPETT-BV421	¥190,000	¥209,000
TB-0174-1	TC-0174-1	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-RLNEVAKNL-PE	¥190,000	¥209,000
TB-0174-2	TC-0174-2	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-RLNEVAKNL-APC	¥190,000	¥209,000
TB-0174-4	TC-0174-4	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-RLNEVAKNL-BV421	¥190,000	¥209,000
TB-0175-1	TC-0175-1	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-KLPDDFTGCV-PE	¥190,000	¥209,000
TB-0175-2	TC-0175-2	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-KLPDDFTGCV-APC	¥190,000	¥209,000
TB-0175-4	TC-0175-4	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-KLPDDFTGCV-BV421	¥190,000	¥209,000
TB-0176-1	TC-0176-1	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-YLQPRTFLL-PE	¥190,000	¥209,000
TB-0176-2	TC-0176-2	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-YLQPRTFLL-APC	¥190,000	¥209,000
TB-0176-4	TC-0176-4	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-YLQPRTFLL-BV421	¥190,000	¥209,000
TB-0177-1	TC-0177-1	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-RLQSLQTYV-PE	¥190,000	¥209,000
TB-0177-2	TC-0177-2	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-RLQSLQTYV-APC	¥190,000	¥209,000
TB-0177-4	TC-0177-4	HLA-A*02:01 SARS-CoV-2 Spike Glycoprotein Tetramer-RLQSLQTYV-BV421	¥190,000	¥209,000
TB-5101-1	TC-5101-1	H-2Db MHV S protein Tetramer-CSLWNGPHL-PE	¥190,000	¥209,000
TB-5101-2	TC-5101-2	H-2Db MHV S protein Tetramer-CSLWNGPHL-APC	¥190,000	¥209,000
TB-5101-4	TC-5101-4	H-2Db MHV S protein Tetramer-CSLWNGPHL-BV421	¥190,000	¥209,000

TB-5102-1	TC-5102-1	H-2Kb Survivin Tetramer-MFFCFKEL-PE	¥190,000	¥209,000
TB-5102-2	TC-5102-2	H-2Kb Survivin Tetramer-MFFCFKEL-APC	¥190,000	¥209,000
TB-5102-4	TC-5102-4	H-2Kb Survivin Tetramer-MFFCFKEL-BV421	¥190,000	¥209,000
TB-5103-1	TC-5103-1	H-2Kb T. gondii tgd057 Tetramer-SVLAFRRL-PE	¥190,000	¥209,000
TB-5103-2	TC-5103-2	H-2Kb T. gondii tgd057 Tetramer-SVLAFRRL-APC	¥190,000	¥209,000
TB-5103-4	TC-5103-4	H-2Kb T. gondii tgd057 Tetramer-SVLAFRRL-BV421	¥190,000	¥209,000
TB-5104-1	TC-5104-1	H-2Kd Tnp03 Tetramer-SYMLQALCI-PE	¥190,000	¥209,000
TB-5104-2	TC-5104-2	H-2Kd Tnp03 Tetramer-SYMLQALCI-APC	¥190,000	¥209,000
TB-5104-4	TC-5104-4	H-2Kd Tnp03 Tetramer-SYMLQALCI-BV421	¥190,000	¥209,000
TB-5105-1	TC-5105-1	H-2Kd MAGE-A3 Tetramer-SYVKVLHHM-PE	¥190,000	¥209,000
TB-5105-2	TC-5105-2	H-2Kd MAGE-A3 Tetramer-SYVKVLHHM-APC	¥190,000	¥209,000
TB-5105-4	TC-5105-4	H-2Kd MAGE-A3 Tetramer-SYVKVLHHM-BV421	¥190,000	¥209,000
TB-5108-1	TC-5108-1	H-2Db Survivin Tetramer-ATFKNWPFL-PE	¥190,000	¥209,000
TB-5108-2	TC-5108-2	H-2Db Survivin Tetramer-ATFKNWPFL-APC	¥190,000	¥209,000
TB-5108-4	TC-5108-4	H-2Db Survivin Tetramer-ATFKNWPFL-BV421	¥190,000	¥209,000
TB-5109-1	TC-5109-1	H-2Db MCMV M45 Tetramer-HGIRNASFI-PE	¥190,000	¥209,000
TB-5109-2	TC-5109-2	H-2Db MCMV M45 Tetramer-HGIRNASFI-APC	¥190,000	¥209,000
TB-5109-4	TC-5109-4	H-2Db MCMV M45 Tetramer-HGIRNASFI-BV421	¥190,000	¥209,000
TB-5110-1	TC-5110-1	H-2Kb HBV HBsAg Tetramer-VWLSVIWM-PE	¥190,000	¥209,000
TB-5110-2	TC-5110-2	H-2Kb HBV HBsAg Tetramer-VWLSVIWM-APC	¥190,000	¥209,000
TB-5110-4	TC-5110-4	H-2Kb HBV HBsAg Tetramer-VWLSVIWM-BV421	¥190,000	¥209,000
TB-5111-1	TC-5111-1	H-2Dd MCMV M164 Tetramer-AGPPRYSRI-PE	¥190,000	¥209,000
TB-5111-2	TC-5111-2	H-2Dd MCMV M164 Tetramer-AGPPRYSRI-APC	¥190,000	¥209,000
TB-5111-4	TC-5111-4	H-2Dd MCMV M164 Tetramer-AGPPRYSRI-BV421	¥190,000	¥209,000
TB-5112-1	TC-5112-1	H-2Ld Mtb Ag85A Tetramer-MPVGQSSF-PE	¥190,000	¥209,000
TB-5112-2	TC-5112-2	H-2Ld Mtb Ag85A Tetramer-MPVGQSSF-APC	¥190,000	¥209,000
TB-5112-4	TC-5112-4	H-2Ld Mtb Ag85A Tetramer-MPVGQSSF-BV421	¥190,000	¥209,000
TB-5113-1	TC-5113-1	H-2Db Adpgk Neopeptide Tetramer-ASMTNMLM-PE	¥190,000	¥209,000
TB-5113-2	TC-5113-2	H-2Db Adpgk Neopeptide Tetramer-ASMTNMLM-APC	¥190,000	¥209,000
TB-5113-4	TC-5113-4	H-2Db Adpgk Neopeptide Tetramer-ASMTNMLM-BV421	¥190,000	¥209,000
TB-5114-1	TC-5114-1	H-2Db REPS1 Neopeptide Tetramer-AQLANDVVL-PE	¥190,000	¥209,000
TB-5114-2	TC-5114-2	H-2Db REPS1 Neopeptide Tetramer-AQLANDVVL-APC	¥190,000	¥209,000
TB-5114-4	TC-5114-4	H-2Db REPS1 Neopeptide Tetramer-AQLANDVVL-BV421	¥190,000	¥209,000
TB-5115-1	TC-5115-1	H-2Db Copine-1 Neopeptide Tetramer-SSPYSLHYL-PE	¥190,000	¥209,000
TB-5115-2	TC-5115-2	H-2Db Copine-1 Neopeptide Tetramer-SSPYSLHYL-APC	¥190,000	¥209,000
TB-5115-4	TC-5115-4	H-2Db Copine-1 Neopeptide Tetramer-SSPYSLHYL-BV421	¥190,000	¥209,000
TB-0008-1	TS-0008-1C	HLA-A*02:01 HIV pol Tetramer-ILKEPVHGV-PE	¥190,000	¥209,000
TB-0008-2	TS-0008-2C	HLA-A*02:01 HIV pol Tetramer-ILKEPVHGV-APC	¥190,000	¥209,000
TB-0008-4	TS-0008-4C	HLA-A*02:01 HIV pol Tetramer-ILKEPVHGV-BV421	¥190,000	¥209,000
TB-0009-1	TS-0009-1C	HLA-A*02:01 Mart-1 Tetramer-ELAGIGILTV-PE	¥190,000	¥209,000
TB-0009-2	TS-0009-2C	HLA-A*02:01 Mart-1 Tetramer-ELAGIGILTV-APC	¥190,000	¥209,000
TB-0009-4	TS-0009-4C	HLA-A*02:01 Mart-1 Tetramer-ELAGIGILTV-BV421	¥190,000	¥209,000
TB-0010-4	TS-0010-4C	HLA-A*02:01 CMV pp65 Tetramer-NLPPMVAIV-BV421	¥190,000	¥209,000



TB-0011-1	TS-0011-1C	HLA-A*02:01 EBV BMLF1 Tetramer-GLCTLVAML-PE	¥190,000	¥209,000
TB-0011-2	TS-0011-2C	HLA-A*02:01 EBV BMLF1 Tetramer-GLCTLVAML-APC	¥190,000	¥209,000
TB-0011-4	TS-0011-4C	HLA-A*02:01 EBV BMLF1 Tetramer-GLCTLVAML-BV421	¥190,000	¥209,000
TB-0012-4	TS-0012-4C	HLA-A*02:01 Influenza-M1 Tetramer-GILGFVFTL-BV421	¥190,000	¥209,000
TB-0013-1	TS-0013-1C	HLA-A*02:01 gp100 (mutant) Tetramer-IMDQVPFSV-PE	¥190,000	¥209,000
TB-0013-2	TS-0013-2C	HLA-A*02:01 gp100 (mutant) Tetramer-IMDQVPFSV-APC	¥190,000	¥209,000
TB-0013-4	TS-0013-4C	HLA-A*02:01 gp100 (mutant) Tetramer-IMDQVPFSV-BV421	¥190,000	¥209,000
TB-0014-1	TS-0014-1C	HLA-A*02:01 gp100 (Wild) Tetramer-ITDQVPFSV-PE	¥190,000	¥209,000
TB-0014-2	TS-0014-2C	HLA-A*02:01 gp100 (Wild) Tetramer-ITDQVPFSV-APC	¥190,000	¥209,000
TB-0014-4	TS-0014-4C	HLA-A*02:01 gp100 (Wild) Tetramer-ITDQVPFSV-BV421	¥190,000	¥209,000
TB-0015-1	TS-0015-1C	HLA-A*02:01 Her-2/neu Tetramer-KIFGSLAFL-PE	¥190,000	¥209,000
TB-0015-2	TS-0015-2C	HLA-A*02:01 Her-2/neu Tetramer-KIFGSLAFL-APC	¥190,000	¥209,000
TB-0015-4	TS-0015-4C	HLA-A*02:01 Her-2/neu E75 Tetramer-KIFGSLAFL-BV421	¥190,000	¥209,000
TB-0016-1	TS-0016-1	HLA-A*02:01 Her-2/neu Tetramer-RLLQETELV-PE	¥190,000	¥209,000
TB-0016-2	TS-0016-2	HLA-A*02:01 Her-2/neu Tetramer-RLLQETELV-APC	¥190,000	¥209,000
TB-0016-4	TS-0016-4	HLA-A*02:01 Her-2/neu Tetramer-RLLQETELV-BV421	¥190,000	¥209,000
TB-0017-1	TS-0017-1	HLA-A*02:01 PR-1 Tetramer-VLQELNVTV-PE	¥190,000	¥209,000
TB-0017-2	TS-0017-2	HLA-A*02:01 PR-1 Tetramer-VLQELNVTV-APC	¥190,000	¥209,000
TB-0017-4	TS-0017-4	HLA-A*02:01 PR-1 Tetramer-VLQELNVTV-BV421	¥190,000	¥209,000
TB-0018-1	TS-0018-1C	HLA-A*02:01 HBV core Tetramer-FLPSDFFPSV-PE	¥190,000	¥209,000
TB-0018-2	TS-0018-2C	HLA-A*02:01 HBV core Tetramer-FLPSDFFPSV-APC	¥190,000	¥209,000
TB-0018-4	TS-0018-4C	HLA-A*02:01 HBV core Tetramer-FLPSDFFPSV-BV421	¥190,000	¥209,000
TB-0019-4	TS-0019-4C	HLA-A*02:01 Tyrosinase Tetramer-YMDGTMSQV-BV421	¥190,000	¥209,000
TB-0020-4	TS-0020-4C	HLA-A*24:02 CMV pp65 Tetramer-QYDPVAALF-BV421	¥190,000	¥209,000
TB-0022-1	TS-0022-1C	HLA-A*24:02 HBV core Tetramer-EYLVSFVGVW-PE	¥190,000	¥209,000
TB-0022-2	TS-0022-2C	HLA-A*24:02 HBV core Tetramer-EYLVSFVGVW-APC	¥190,000	¥209,000
TB-0022-4	TS-0022-4C	HLA-A*24:02 HBV core Tetramer-EYLVSFVGVW-BV421	¥190,000	¥209,000
TB-0023-4	TS-0023-4C	HLA-A*24:02 HBV pol Tetramer-KYTSFPWLL-BV421	¥190,000	¥209,000
TB-0024-4	TS-0024-4C	HLA-A*01:01 CMV pp50 Tetramer-VTEHDTLLY-BV421	¥190,000	¥209,000
TB-0025-1	TS-0025-1C	HLA-B*07:02 CMV pp65 Tetramer-TPRVTGGGAM-PE	¥190,000	¥209,000
TB-0025-2	TS-0025-2C	HLA-B*07:02 CMV pp65 Tetramer-TPRVTGGGAM-APC	¥190,000	¥209,000
TB-0025-4	TS-0025-4C	HLA-B*07:02 CMV pp65 Tetramer-TPRVTGGGAM-BV421	¥190,000	¥209,000
TB-0026-4	TS-0026-4C	HLA-B*08:01 CMV IE1 Tetramer-ELRRKMMYM-BV421	¥190,000	¥209,000
TB-0027-1	TS-0027-1C	HLA-B*35:01 CMV pp65 Tetramer-IPSINVHHY-PE	¥190,000	¥209,000
TB-0027-2	TS-0027-2C	HLA-B*35:01 CMV pp65 Tetramer-IPSINVHHY-APC	¥190,000	¥209,000
TB-0027-4	TS-0027-4C	HLA-B*35:01 CMV pp65 Tetramer-IPSINVHHY-BV421	¥190,000	¥209,000
TB-0029-1	TS-0029-1C	HLA-A*02:01 Negative Tetramer-PE	¥190,000	¥209,000
TB-0029-2	TS-0029-2C	HLA-A*02:01 Negative Tetramer-APC	¥190,000	¥209,000
TB-0029-4	TS-0029-4C	HLA-A*02:01 Negative Tetramer-BV421	¥190,000	¥209,000
TB-0031-1	TS-0031-1	HLA-A*02:01 HPV16 E7 Tetramer-YMLDLQPET-PE	¥190,000	¥209,000
TB-0031-2	TS-0031-2	HLA-A*02:01 HPV16 E7 Tetramer-YMLDLQPET-APC	¥190,000	¥209,000
TB-0031-4	TS-0031-4	HLA-A*02:01 HPV16 E7 Tetramer-YMLDLQPET-BV421	¥190,000	¥209,000
TB-0035-1	TS-0035-1C	HLA-A*02:01 gp100 154-162 Tetramer-KTWGQYVQV-PE	¥190,000	¥209,000



TB-0035-2	TS-0035-2C	HLA-A*02:01 gp100 154-162 Tetramer-KTWGQYWQV-APC	¥190,000	¥209,000
TB-0035-4	TS-0035-4C	HLA-A*02:01 gp100 154-162 Tetramer-KTWGQYWQV-BV421	¥190,000	¥209,000
TB-5001-4	TS-5001-4C	H-2Kb OVA Tetramer-SIINFEKL-BV421	¥190,000	¥209,000
TB-5002-1	TS-5002-1C	H-2Db LCMV Tetramer-KAVYNFATC-PE	¥190,000	¥209,000
TB-5002-4	TS-5002-4C	H-2Db LCMV gp33 Tetramer-KAVYNFATC-BV421	¥190,000	¥209,000
TB-5004-1	TS-5004-1C	H-2Kb TRP-2 Tetramer-SVYDFVWL-PE	¥190,000	¥209,000
TB-5004-2	TS-5004-2C	H-2Kb TRP-2 Tetramer-SVYDFVWL-APC	¥190,000	¥209,000
TB-5004-4	TS-5004-4C	H-2Kb TRP-2 Tetramer-SVYDFVWL-BV421	¥190,000	¥209,000
TB-5007-1	TS-5007-1	H-2Kd HIV gag Tetramer-AMQMLKETI-PE	¥190,000	¥209,000
TB-5007-2	TS-5007-2	H-2Kd HIV gag Tetramer-AMQMLKETI-APC	¥190,000	¥209,000
TB-5007-4	TS-5007-4	H-2Kd HIV gag Tetramer-AMQMLKETI-BV421	¥190,000	¥209,000
TB-5008-1	TS-5008-1C	H-2Db HPV16 E7 Tetramer-RAHYNIVTF-PE	¥190,000	¥209,000
TB-5008-2	TS-5008-2C	H-2Db HPV16 E7 Tetramer-RAHYNIVTF-APC	¥190,000	¥209,000
TB-5008-4	TS-5008-4C	H-2Db HPV16 E7 Tetramer-RAHYNIVTF-BV421	¥190,000	¥209,000
TB-5009-1	TS-5009-1	H-2Db LCMV gp276-286 Tetramer-SGVENPGGYCL-PE	¥190,000	¥209,000
TB-5009-2	TS-5009-2	H-2Db LCMV gp276-286 Tetramer-SGVENPGGYCL-APC	¥190,000	¥209,000
TB-5009-4	TS-5009-4	H-2Db LCMV gp276-286 Tetramer-SGVENPGGYCL-BV421	¥190,000	¥209,000
TB-5010-1	TS-5010-1C	H-2Kb LCMV gp34-41 Tetramer-AVYNFATC-PE	¥190,000	¥209,000
TB-5010-2	TS-5010-2C	H-2Kb LCMV gp34-41 Tetramer-AVYNFATC-APC	¥190,000	¥209,000
TB-5010-4	TS-5010-4C	H-2Kb LCMV gp34-41 Tetramer-AVYNFATC-BV421	¥190,000	¥209,000
TB-5011-1	TS-5011-1	H-2Kb LCMV gp34-43 Tetramer-AVYNFATCGI-PE	¥190,000	¥209,000
TB-5011-2	TS-5011-2	H-2Kb LCMV gp34-43 Tetramer-AVYNFATCGI-APC	¥190,000	¥209,000
TB-5011-4	TS-5011-4	H-2Kb LCMV gp34-43 Tetramer-AVYNFATCGI-BV421	¥190,000	¥209,000
TB-5012-1	TS-5012-1	H-2Kb LCMV gp118-125 Tetramer-ISHNFCNL-PE	¥190,000	¥209,000
TB-5012-2	TS-5012-2	H-2Kb LCMV gp118-125 Tetramer-ISHNFCNL-APC	¥190,000	¥209,000
TB-5012-4	TS-5012-4	H-2Kb LCMV gp118-125 Tetramer-ISHNFCNL-BV421	¥190,000	¥209,000
TB-5014-1	TS-5014-1	H-2Kb LCMV L protein Tetramer-LEYDFNKL-PE	¥190,000	¥209,000
TB-5014-2	TS-5014-2	H-2Kb LCMV L protein Tetramer-LEYDFNKL-APC	¥190,000	¥209,000
TB-5014-4	TS-5014-4	H-2Kb LCMV L protein Tetramer-LEYDFNKL-BV421	¥190,000	¥209,000
TB-5015-1	TS-5015-1	H-2Kb LCMV NP205-212 Tetramer-YTVKYPNL-PE	¥190,000	¥209,000
TB-5015-2	TS-5015-2	H-2Kb LCMV NP205-212 Tetramer-YTVKYPNL-APC	¥190,000	¥209,000
TB-5015-4	TS-5015-4	H-2Kb LCMV NP205-212 Tetramer-YTVKYPNL-BV421	¥190,000	¥209,000
TB-5016-1	TS-5016-1	H-2Db MoMSV Tetramer-(Abu)(Abu)L(Abu)LTVFL-PE	¥190,000	¥209,000
TB-5016-2	TS-5016-2	H-2Db MoMSV Tetramer-(Abu)(Abu)L(Abu)LTVFL-APC	¥190,000	¥209,000
TB-5016-4	TS-5016-4	H-2Db MoMSV Tetramer-(Abu)(Abu)L(Abu)LTVFL-BV421	¥190,000	¥209,000
TB-5017-1	TS-5017-1C	H-2Db SIV gag Tetramer-AAVKNWMTQTL-PE	¥190,000	¥209,000
TB-5017-2	TS-5017-2C	H-2Db SIV gag Tetramer-AAVKNWMTQTL-APC	¥190,000	¥209,000
TB-5017-4	TS-5017-4C	H-2Db SIV gag Tetramer-AAVKNWMTQTL-BV421	¥190,000	¥209,000
TB-5018-1	TS-5018-1C	H-2Db RSV M187-195 Tetramer-NAITNAKII-PE	¥190,000	¥209,000
TB-5018-2	TS-5018-2C	H-2Db RSV M187-195 Tetramer-NAITNAKII-APC	¥190,000	¥209,000
TB-5018-4	TS-5018-4C	H-2Db RSV M187-195 Tetramer-NAITNAKII-BV421	¥190,000	¥209,000
TB-7300-K1	TS-7300-1F	QuickSwitch Quant HLA-A*02:01 Tetramer Kit-PE (frozen)	¥157,000	¥188,000
	TS-7300-1R	QuickSwitch Quant HLA-A*02:01 Tetramer Kit-PE (2-8°C)		

TB-7300-K2	TS-7300-2F	QuickSwitch Quant HLA-A*02:01 Tetramer Kit-APC (frozen)	¥157,000	¥188,000
	TS-7300-2R	QuickSwitch Quant HLA-A*02:01 Tetramer Kit-APC (2-8°C)		
TB-7300-KM	TS-7300-SF	QuickSwitch Quant HLA-A*02:01 Peptide Screening Kit (frozen)	¥330,000	¥330,000 (変更なし)
	TS-7300-SR	QuickSwitch Quant HLA-A*02:01 Peptide Screening Kit (2-8°C)		
TB-7302-K1	TS-7302-1F	QuickSwitch Quant HLA-A*24:02 Tetramer Kit-PE (frozen)	¥157,000	¥188,000
	TS-7302-1R	QuickSwitch Quant HLA-A*24:02 Tetramer Kit-PE (2-8°C)		
TB-7302-K4	TS-7302-4F	QuickSwitch Quant HLA-A*24:02 Tetramer Kit-BV421 (frozen)	¥157,000	¥188,000
	TS-7302-4R	QuickSwitch Quant HLA-A*24:02 Tetramer Kit-BV421 (2-8°C)		
TB-7304-K1	TS-7304-1F	QuickSwitch Quant HLA-A*11:01 Tetramer Kit-PE (frozen)	¥157,000	¥188,000
	TS-7304-1R	QuickSwitch Quant HLA-A*11:01 Tetramer Kit-PE (2-8°C)		
TB-7304-K2	TS-7304-2F	QuickSwitch Quant HLA-A*11:01 Tetramer Kit-APC (frozen)	¥157,000	¥188,000
	TS-7304-2R	QuickSwitch Quant HLA-A*11:01 Tetramer Kit-APC (2-8°C)		
TB-7304-K4	TS-7304-4F	QuickSwitch Quant HLA-A*11:01 Tetramer Kit-BV421 (frozen)	¥157,000	¥188,000
	TS-7304-4R	QuickSwitch Quant HLA-A*11:01 Tetramer Kit-BV421 (2-8°C)		
TB-7306-K1	TS-7306-1F	QuickSwitch Quant HLA-A*03:01 Tetramer Kit-PE (frozen)	¥157,000	¥188,000
	TS-7306-1R	QuickSwitch Quant HLA-A*03:01 Tetramer Kit-PE (2-8°C)		
TB-7306-K2	TS-7306-2F	QuickSwitch Quant HLA-A*03:01 Tetramer Kit-APC (frozen)	¥157,000	¥188,000
	TS-7306-2R	QuickSwitch Quant HLA-A*03:01 Tetramer Kit-APC (2-8°C)		
TB-7306-K4	TS-7306-4F	QuickSwitch Quant HLA-A*03:01 Tetramer Kit-BV421 (frozen)	¥157,000	¥188,000
	TS-7306-4R	QuickSwitch Quant HLA-A*03:01 Tetramer Kit-BV421 (2-8°C)		
TB-7308-K1	TS-7308-1F	QuickSwitch Quant HLA-A*01:01 Tetramer Kit-PE (frozen)	¥157,000	¥188,000
	TS-7308-1R	QuickSwitch Quant HLA-A*01:01 Tetramer Kit-PE (2-8°C)		
TB-7308-K2	TS-7308-2F	QuickSwitch Quant HLA-A*01:01 Tetramer Kit-APC (frozen)	¥157,000	¥188,000
	TS-7308-2R	QuickSwitch Quant HLA-A*01:01 Tetramer Kit-APC (2-8°C)		
TB-7308-K4	TS-7308-4F	QuickSwitch Quant HLA-A*01:01 Tetramer Kit-BV421 (frozen)	¥157,000	¥188,000
	TS-7308-4R	QuickSwitch Quant HLA-A*01:01 Tetramer Kit-BV421 (2-8°C)		
TB-7310-K1	TS-7310-1F	QuickSwitch Quant HLA-B*07:02 Tetramer Kit-PE (frozen)	¥157,000	¥188,000
	TS-7310-1R	QuickSwitch Quant HLA-B*07:02 Tetramer Kit-PE (2-8°C)		
TB-7310-K2	TS-7310-2F	QuickSwitch Quant HLA-B*07:02 Tetramer Kit-APC (frozen)	¥157,000	¥188,000
	TS-7310-2R	QuickSwitch Quant HLA-B*07:02 Tetramer Kit-APC (2-8°C)		
TB-7310-K4	TS-7310-4F	QuickSwitch Quant HLA-B*07:02 Tetramer Kit-BV421 (frozen)	¥157,000	¥188,000
	TS-7310-4R	QuickSwitch Quant HLA-B*07:02 Tetramer Kit-BV421 (2-8°C)		
TB-7400-K1	TS-7400-1F	QuickSwitch Quant H-2Kb Tetramer Kit-PE (frozen)	¥157,000	¥188,000
	TS-7400-1R	QuickSwitch Quant H-2Kb Tetramer Kit-PE (2-8°C)		
TB-7400-K2	TS-7400-2F	QuickSwitch Quant H-2Kb Tetramer Kit-APC (frozen)	¥157,000	¥188,000
	TS-7400-2R	QuickSwitch Quant H-2Kb Tetramer Kit-APC (2-8°C)		
TB-7400-K4	TS-7400-4F	QuickSwitch Quant H-2Kb Tetramer Kit-BV421 (frozen)	¥157,000	¥188,000
	TS-7400-4R	QuickSwitch Quant H-2Kb Tetramer Kit-BV421 (2-8°C)		
TB-M011-1	TS-M011-1	HLA-A*02:01 NY-ESO-1 Tetramer-SLLMWITQC-PE	¥190,000	¥209,000
TB-M011-2	TS-M011-2	HLA-A*02:01 NY-ESO-1 Tetramer-SLLMWITQC-APC	¥190,000	¥209,000
TB-M011-4	TS-M011-4	HLA-A*02:01 NY-ESO-1 Tetramer-SLLMWITQC-BV421	¥190,000	¥209,000
TB-M014-4	TS-M014-4	HLA-A*24:02 modified WT1 Tetramer-CYTWNQMNL-BV421	¥190,000	¥209,000
TB-M016-4	TS-M016-4	HLA-A*02:01 WT-1126-134 Tetramer-RMFNPAPYL-BV421	¥190,000	¥209,000

TB-M017-4	TS-M017-4	HLA-A*02:01 HTLV-1 Tax11-19 Tetramer-LLFGYPVYV-BV421	¥190,000	¥209,000
TB-M027-4	TS-M027-4	HLA-A*02:01 HIV gag Tetramer-SLYNTVATL-BV421	¥190,000	¥209,000
TB-M037-1	TS-M037-1	HLA-B*35:01 EBV BZLF1 54-64 Tetramer-EPLPQQQLTAY-PE	¥190,000	¥209,000
TB-M037-2	TS-M037-2	HLA-B*35:01 EBV BZLF1 54-64 Tetramer-EPLPQQQLTAY-APC	¥190,000	¥209,000
TB-M038-1	TS-M038-1	HLA-B*35:01 EBV LMP2 1-9 Tetramer-MGSLEMVPM-PE	¥190,000	¥209,000
TB-M038-2	TS-M038-2	HLA-B*35:01 EBV LMP2 1-9 Tetramer-MGSLEMVPM-APC	¥190,000	¥209,000
TB-M039-1	TS-M039-1	HLA-A*02:01 HCV NS3 1073-1081 Tetramer-CINGVCWTV-PE	¥190,000	¥209,000
TB-M039-2	TS-M039-2	HLA-A*02:01 HCV NS3 1073-1081 Tetramer-CINGVCWTV-APC	¥190,000	¥209,000
TB-M039-4	TS-M039-4	HLA-A*02:01 HCV NS3 1073-1081 Tetramer-CINGVCWTV-BV421	¥190,000	¥209,000
TB-M040-1	TS-M040-1	HLA-A*02:01 HCV NS3 1406-1415 Tetramer-KLVALGINAV-PE	¥190,000	¥209,000
TB-M040-2	TS-M040-2	HLA-A*02:01 HCV NS3 1406-1415 Tetramer-KLVALGINAV-APC	¥190,000	¥209,000
TB-M046-1	TS-M046-1	HLA-B*35:01 Influenza NP Tetramer-LPFEKSTVM-PE	¥190,000	¥209,000
TB-M046-2	TS-M046-2	HLA-B*35:01 Influenza NP Tetramer-LPFEKSTVM-APC	¥190,000	¥209,000
TB-M047-1	TS-M047-1	HLA-A*02:01 HPV16 E6 Tetramer-KLPQLCTEL-PE	¥190,000	¥209,000
TB-M047-2	TS-M047-2	HLA-A*02:01 HPV16 E6 Tetramer-KLPQLCTEL-APC	¥190,000	¥209,000
TB-M047-4	TS-M047-4	HLA-A*02:01 HPV16 E6 Tetramer-KLPQLCTEL-BV421	¥190,000	¥209,000
TB-M048-4	TS-M048-4	HLA-A*02:01 HPV16 E7 Tetramer-YMLDLQPETT-BV421	¥190,000	¥209,000
TB-M051-1	TS-M051-1	HLA-A*02:01 HBV env 335-343 Tetramer-WLSLLVPFV-PE	¥190,000	¥209,000
TB-M051-2	TS-M051-2	HLA-A*02:01 HBV env 335-343 Tetramer-WLSLLVPFV-APC	¥190,000	¥209,000
TB-M053-1	TS-M053-1	HLA-A*02:01 HBV pol Tetramer-FLLSLGIHL-PE	¥190,000	¥209,000
TB-M053-2	TS-M053-2	HLA-A*02:01 HBV pol Tetramer-FLLSLGIHL-APC	¥190,000	¥209,000
TB-M055-1	TS-M055-1	HLA-B*35:01 HIV RT Tetramer-NPDIVYQY-PE	¥190,000	¥209,000
TB-M057-1	TS-M057-1	HLA-A*02:01 CMV IE1 316-324 Tetramer-VLEETSVML-PE	¥190,000	¥209,000
TB-M057-2	TS-M057-2	HLA-A*02:01 CMV IE1 316-324 Tetramer-VLEETSVML-APC	¥190,000	¥209,000
TB-M059-1	TS-M059-1	HLA-A*02:01 Adenovirus 11 Hexon914-922 Tetramer-LLFEVFDVV-PE	¥190,000	¥209,000
TB-M059-2	TS-M059-2	HLA-A*02:01 Adenovirus 11 Hexon914-922 Tetramer-LLFEVFDVV-APC	¥190,000	¥209,000
TB-M067-1	TS-M067-1	HLA-B*35:01 AdV Hexon 320-329 Tetramer-MPNRPNYIAF-PE	¥190,000	¥209,000
TB-M067-2	TS-M067-2	HLA-B*35:01 AdV Hexon 320-329 Tetramer-MPNRPNYIAF-APC	¥190,000	¥209,000
TB-M068-1	TS-M068-1	HLA-B*35:01 AdV Hexon 705-713 Tetramer-IPYLDGTFY-PE	¥190,000	¥209,000
TB-M068-2	TS-M068-2	HLA-B*35:01 AdV Hexon 705-713 Tetramer-IPYLDGTFY-APC	¥190,000	¥209,000
TB-M070-1	TS-M070-1	HLA-A*02:01 MAGE-A1 Tetramer-KVLEYVIKV-PE	¥190,000	¥209,000
TB-M070-2	TS-M070-2	HLA-A*02:01 MAGE-A1 Tetramer-KVLEYVIKV-APC	¥190,000	¥209,000
TB-M070-4	TS-M070-4	HLA-A*02:01 MAGE-A1 Tetramer-KVLEYVIKV-BV421	¥190,000	¥209,000
TB-M072-1	TS-M072-1	HLA-A*02:01 MAGE-A2 Tetramer-YLQLVFGIEV-PE	¥190,000	¥209,000
TB-M072-2	TS-M072-2	HLA-A*02:01 MAGE-A2 Tetramer-YLQLVFGIEV-APC	¥190,000	¥209,000
TB-M074-1	TS-M074-1	HLA-A*01:01 MAGE-A3 Tetramer-EVDPIGHLY-PE	¥190,000	¥209,000
TB-M074-2	TS-M074-2	HLA-A*01:01 MAGE-A3 Tetramer-EVDPIGHLY-APC	¥190,000	¥209,000
TB-M074-4	TS-M074-4	HLA-A*01:01 MAGE-A3 Tetramer-EVDPIGHLY-BV421	¥190,000	¥209,000
TB-M075-1	TS-M075-1	HLA-A*02:01 MAGE-A3 112-120 Tetramer-KVAELVHFL-PE	¥190,000	¥209,000
TB-M075-2	TS-M075-2	HLA-A*02:01 MAGE-A3 112-120 Tetramer-KVAELVHFL-APC	¥190,000	¥209,000
TB-M076-1	TS-M076-1	HLA-A*02:01 MAGE-A3 271-279 Tetramer-FLWGPRALV-PE	¥190,000	¥209,000
TB-M076-2	TS-M076-2	HLA-A*02:01 MAGE-A3 271-279 Tetramer-FLWGPRALV-APC	¥190,000	¥209,000
TB-M082-1	TS-M082-1	HLA-A*02:01 gp100 Tetramer-YLEPGPVTA-PE	¥190,000	¥209,000

TB-M082-2	TS-M082-2	HLA-A*02:01 gp100 Tetramer-YLEPGPVTA-APC	¥190,000	¥209,000
TB-M082-4	TS-M082-4	HLA-A*02:01 gp100 Tetramer-YLEPGPVTA-BV421	¥190,000	¥209,000
TB-M088-1	TS-M088-1	HLA-A*02:01 MUC1 Tetramer-LLLLTVLTV-PE	¥190,000	¥209,000
TB-M088-2	TS-M088-2	HLA-A*02:01 MUC1 Tetramer-LLLLTVLTV-APC	¥190,000	¥209,000
TB-M099-1	TS-M099-1	HLA-B*07:02 CMV pp65 Tetramer-RPHERNGFTVL-PE	¥190,000	¥209,000
TB-M099-2	TS-M099-2	HLA-B*07:02 CMV pp65 Tetramer-RPHERNGFTVL-APC	¥190,000	¥209,000
TB-M105-1	TS-M105-1	HLA-A*02:01 NY-ESO-1 C9V Tetramer-SLLMWITQV-PE	¥190,000	¥209,000
TB-M105-2	TS-M105-2	HLA-A*02:01 NY-ESO-1 C9V Tetramer-SLLMWITQV-APC	¥190,000	¥209,000
TB-M105-4	TS-M105-4	HLA-A*02:01 NY-ESO-1 C9V Tetramer-SLLMWITQV-BV421	¥190,000	¥209,000
TB-M106-1	TS-M106-1	HLA-B*35:01 HIV nef 74-81 Tetramer-VPLRPMTY-PE	¥190,000	¥209,000
TB-M106-2	TS-M106-2	HLA-B*35:01 HIV nef 74-81 Tetramer-VPLRPMTY-APC	¥190,000	¥209,000
TB-M107-1	TS-M107-1	HLA-A*02:01 PAP 299-307 Tetramer-ALDVYNGLL-PE	¥190,000	¥209,000
TB-M107-2	TS-M107-2	HLA-A*02:01 PAP 299-307 Tetramer-ALDVYNGLL-APC	¥190,000	¥209,000
TB-M120-1	TS-M120-1	HLA-A*02:01 PSA141-150 Tetramer-FLTPKKLQCV-PE	¥190,000	¥209,000
TB-M120-2	TS-M120-2	HLA-A*02:01 PSA141-150 Tetramer-FLTPKKLQCV-APC	¥190,000	¥209,000
TB-M122-1	TS-M122-1	HLA-A*02:01 VZV IE62 Tetramer-ALWALPHAA-PE	¥190,000	¥209,000
TB-M122-2	TS-M122-2	HLA-A*02:01 VZV IE62 Tetramer-ALWALPHAA-APC	¥190,000	¥209,000
TB-M505-4	TS-M505-4	H-2Db human gp 100 Tetramer-KVPRNQDWL-BV421	¥190,000	¥209,000
TB-M506-4	TS-M506-4	H-2Kd RSV M2 Tetramer-SYIGSINNI-BV421	¥190,000	¥209,000
TB-M508-4	TS-M508-4	H-2Db Influenza NP Tetramer-ASNENMETM-BV421	¥190,000	¥209,000
TB-M513-4	TS-M513-4	H-2Db LCMV NP369 Tetramer-FQPQNGQFI-BV421	¥190,000	¥209,000
TB-M524-2	TS-M524-2	H-2Db HY Uty Tetramer-WMHHNMDLI-APC	¥190,000	¥209,000
TB-M524-4	TS-M524-4	H-2Db HY Uty Tetramer-WMHHNMDLI-BV421	¥190,000	¥209,000
TB-M534-1	TS-M534-1	H-2Kd Influenza NP Tetramer-TYQRTRALV-PE	¥190,000	¥209,000
TB-M534-2	TS-M534-2	H-2Kd Influenza NP Tetramer-TYQRTRALV-APC	¥190,000	¥209,000
TB-M536-1	TS-M536-1	H-2Dd HIV env Tetramer-IGPGRAFYA-PE	¥190,000	¥209,000
TB-M536-2	TS-M536-2	H-2Dd HIV env Tetramer-IGPGRAFYA-APC	¥190,000	¥209,000
TB-M537-1	TS-M537-1	H-2Kb HBV core Tetramer-MGLKFRQL-PE	¥190,000	¥209,000
TB-M537-2	TS-M537-2	H-2Kb HBV core Tetramer-MGLKFRQL-APC	¥190,000	¥209,000
TB-M538-1	TS-M538-1	H-2Kb VACV B8R Tetramer-TSYKFESV-PE	¥190,000	¥209,000
TB-M538-2	TS-M538-2	H-2Kb VACV B8R Tetramer-TSYKFESV-APC	¥190,000	¥209,000
TB-M539-1	TS-M539-1	H-2Db SV40 large T Ag 206-215 Tetramer-SAINNYAQKL-PE	¥190,000	¥209,000
TB-M539-2	TS-M539-2	H-2Db SV40 large T Ag 206-215 Tetramer-SAINNYAQKL-APC	¥190,000	¥209,000
TB-M539-4	TS-M539-4	H-2Db SV40 large T Ag 206-215 Tetramer-SAINNYAQKL-BV421	¥190,000	¥209,000
TB-M547-1	TS-M547-1	H-2Kd Plasmodium CSP Tetramer-SYVPSAEQI-PE	¥190,000	¥209,000
TB-M547-2	TS-M547-2	H-2Kd Plasmodium CSP Tetramer-SYVPSAEQI-APC	¥190,000	¥209,000
TB-M552-1	TS-M552-1	H-2Kd IGRP Tetramer-VYLKTNVFL-PE	¥190,000	¥209,000
TB-M552-2	TS-M552-2	H-2Kd IGRP Tetramer-VYLKTNVFL-APC	¥190,000	¥209,000
TB-M553-1	TS-M553-1	H-2Kd NRP-V7 Tetramer-KYNKANVFL-PE	¥190,000	¥209,000
TB-M553-2	TS-M553-2	H-2Kd NRP-V7 Tetramer-KYNKANVFL-APC	¥190,000	¥209,000
TB-M562-1	TS-M562-1	H-2Kb mTERT Tetramer-VGRNFTNL-PE	¥190,000	¥209,000
TB-M562-2	TS-M562-2	H-2Kb mTERT Tetramer-VGRNFTNL-APC	¥190,000	¥209,000

## Appendix 2 : New products

Code No.	Product Name	Size	Price (w/o tax)
TC-0117-1	HLA-A*24:02 EBV LMP-2 Tetramer-TYGPVFMCL-PE	50 tests	¥209,000
TC-0117-2	HLA-A*24:02 EBV LMP-2 Tetramer-TYGPVFMCL-APC	50 tests	¥209,000
TC-0165-1	HLA-A*02:01 Influenza M1 Tetramer-ILGFVFTLTV-PE	50 tests	¥209,000
TC-0165-2	HLA-A*02:01 Influenza M1 Tetramer-ILGFVFTLTV-APC	50 tests	¥209,000
TC-5116-1	H-2Kb RPL18 Neoepitope Tetramer-KILTFDRL-PE	50 tests	¥209,000
TC-5116-2	H-2Kb RPL18 Neoepitope Tetramer-KILTFDRL-APC	50 tests	¥209,000
TC-5116-4	H-2Kb RPL18 Neoepitope Tetramer-KILTFDRL-BV421	50 tests	¥209,000
TS-5002-2C	H-2Db LCMV Tetramer-KAVYNFATC-APC	50 tests	¥209,000
TS-7300-4F	QuickSwitch Quant HLA-A*02:01 Tetramer Kit-BV421 (frozen)	25 µg	¥188,000
TS-7300-4R	QuickSwitch Quant HLA-A*02:01 Tetramer Kit-BV421 (2-8°C)	25 µg	
TS-7300-MF	QuickSwitch Quant HLA-A*02:01 Monomer Kit (frozen)	25 µg	¥188,000
TS-7300-MR	QuickSwitch Quant HLA-A*02:01 Monomer Kit (2-8°C)	25 µg	
TS-7302-2F	QuickSwitch Quant HLA-A*24:02 Tetramer Kit-APC (frozen)	25 µg	¥188,000
TS-7302-2R	QuickSwitch Quant HLA-A*24:02 Tetramer Kit-APC (2-8°C)	25 µg	
TS-M007-4	HLA-A*24:02 Negative (HIV env) Tetramer-RYLRDQQLL-BV421	50 tests	¥209,000
TS-M008-4	H-2Kb Negative Tetramer-SIYRYGL-BV421	50 tests	¥209,000
TS-M032-4	HLA-A*02:01 EBV LMP2 Tetramer-CLGGLTMMV-BV421	50 tests	¥209,000
TS-M033-4	HLA-A*03:01 EBV EBNA3A 603-611 Tetramer-RLRAEAQVK-BV421	50 tests	¥209,000
TS-M037-4	HLA-B*35:01 EBV BZLF1 54-64 Tetramer-EPLPQGQLTAY-BV421	50 tests	¥209,000
TS-M038-4	HLA-B*35:01 EBV LMP2 1-9 Tetramer-MGSLEMPM-BV421	50 tests	¥209,000
TS-M040-4	HLA-A*02:01 HCV NS3 1406-1415 Tetramer-KLVALGINAV-BV421	50 tests	¥209,000
TS-M043-4	HLA-A*02:01 HCV NS5B 2594-2602 Tetramer-ALYDVVSKL-BV421	50 tests	¥209,000
TS-M046-4	HLA-B*35:01 Influenza NP Tetramer-LPFEKSTVM-BV421	50 tests	¥209,000
TS-M051-4	HLA-A*02:01 HBV env 335-343 Tetramer-WLSLLVPFV-BV421	50 tests	¥209,000
TS-M053-4	HLA-A*02:01 HBV pol Tetramer-FLLSLGIHL-BV421	50 tests	¥209,000
TS-M055-2	HLA-B*35:01 HIV RT Tetramer-NPDIVYQY-APC	50 tests	¥209,000
TS-M055-4	HLA-B*35:01 HIV RT Tetramer-NPDIVYQY-BV421	50 tests	¥209,000
TS-M057-4	HLA-A*02:01 CMV IE1 316-324 Tetramer-VLEETSVML-BV421	50 tests	¥209,000
TS-M059-4	HLA-A*02:01 AdV 11 Hexon 914-922 Tetramer-LLFEVFDVV-BV421	50 tests	¥209,000
TS-M067-4	HLA-B*35:01 AdV Hexon 320-329 Tetramer-MPNRPNYIAF-BV421	50 tests	¥209,000
TS-M068-4	HLA-B*35:01 AdV Hexon 705-713 Tetramer-IPYLDGTFY-BV421	50 tests	¥209,000
TS-M072-4	HLA-A*02:01 MAGE-A2 Tetramer-YLQLVFGIEV-BV421	50 tests	¥209,000
TS-M075-4	HLA-A*02:01 MAGE-A3 112-120 Tetramer-KVAELVHFL-BV421	50 tests	¥209,000
TS-M076-4	HLA-A*02:01 MAGE-A3 271-279 Tetramer-FLWGPRLV-BV421	50 tests	¥209,000
TS-M078-4	HLA-A*02:01 MAGE-A10 Tetramer-GLYDGMEHL-BV421	50 tests	¥209,000
TS-M088-4	HLA-A*02:01 MUC1 Tetramer-LLLLTVLTV-BV421	50 tests	¥209,000
TS-M099-4	HLA-B*07:02 CMV pp65 Tetramer-RPHERNGFTVL-BV421	50 tests	¥209,000

TS-M106-4	HLA-B*35:01 HIV nef 74-81 Tetramer-VPLRPMTY-BV421	50 tests	¥209,000
TS-M107-4	HLA-A*02:01 PAP 299-307 Tetramer-ALDVYNGLL-BV421	50 tests	¥209,000
TS-M115-4	HLA-A*02:01 hTERT Tetramer-ILAKFLHWL-BV421	50 tests	¥209,000
TS-M120-4	HLA-A*02:01 PSA 141-150 Tetramer-FLTPKKLQCV-BV421	50 tests	¥209,000
TS-M122-4	HLA-A*02:01 VZV IE62 Tetramer-ALWALPHAA-BV421	50 tests	¥209,000
TS-M124-4	HLA-A*03:01 EBV BRLF1 Tetramer-RVRAYTYSK-BV421	50 tests	¥209,000
TS-M142-4	HLA-B*07:02 EBV EBNA3A Tetramer-RPPIFIRRL-BV421	50 tests	¥209,000
TS-M507-4	H-2Kb MuLV p15E Tetramer-KSPWFRTL-BV421	50 tests	¥209,000
TS-M511-4	H-2Ld $\beta$ -galactosidase Tetramer-TPHPARIGL-BV421	50 tests	¥209,000
TS-M512-4	H-2Db LCMV gp33 (C9M) Tetramer-KAVYNFATM-BV421	50 tests	¥209,000
TS-M514-2	H-2Ld LCMV NP 118-126 Tetramer-RPQASGVYM-APC	50 tests	¥209,000
TS-M514-4	H-2Ld LCMV NP 118-126 Tetramer-RPQASGVYM-BV421	50 tests	¥209,000
TS-M520-4	H-2Kd Influenza HA Tetramer-IYSTVASSL-BV421	50 tests	¥209,000
TS-M521-4	H-2Ld MuLV gp70 Tetramer-SPSYVYHQF-BV421	50 tests	¥209,000
TS-M528-4	H-2Db Influenza PA Tetramer-SSLENFRAYV-BV421	50 tests	¥209,000
TS-M534-4	H-2Kd Influenza NP Tetramer-TYQRTRALV-BV421	50 tests	¥209,000
TS-M536-4	H-2Dd HIV env Tetramer-IGPGRAFYA-BV421	50 tests	¥209,000
TS-M537-4	H-2Kb HBV core Tetramer-MGLKFRQL-BV421	50 tests	¥209,000
TS-M538-4	H-2Kb VACV B8R Tetramer-TSYKFESV-BV421	50 tests	¥209,000
TS-M547-4	H-2Kd Plasmodium CSP Tetramer-SYVPSAEQI-BV421	50 tests	¥209,000
TS-M552-4	H-2Kd IGRP Tetramer-VYLKTNVFL-BV421	50 tests	¥209,000
TS-M553-4	H-2Kd NRP-V7 Tetramer-KYNKANVFL-BV421	50 tests	¥209,000
TS-M561-4	H-2Db hPSA Tetramer-HCIRNKSVIL-BV421	50 tests	¥209,000
TS-M562-4	H-2Kb mTERT Tetramer-VGRNFTNL-BV421	50 tests	¥209,000
TS-M801-4	HLA-DRB1*01:01 CLIP 103-117 Tetramer-PVSKMRMATPLLMQA-BV421	20 tests	¥198,000
TS-M802-4	HLA-DRB1*01:01 HIV gag 295-307 Tetramer-DYVDRFYKTLRAE-BV421	20 tests	¥198,000
TS-M804-4	HLA-DRB1*01:01 Influenza HA 306-318 Tetramer-PKYVKQNTLKLAT-BV421	20 tests	¥198,000
TS-M809-4	HLA-DRB1*04:01 CLIP 103-117 Tetramer-PVSKMRMATPLLMQA-BV421	20 tests	¥198,000
TS-M810-4	HLA-DRB1*04:01 Influenza HA 306-318 Tetramer-PKYVKQNTLKLAT-BV421	20 tests	¥198,000
TS-ME03-2	HLA-E*01:01 HLA-A leader 3-11 Tetramer-VMAPRTLVL-APC	50 tests	¥209,000
TS-ME03-4	HLA-E*01:01 HLA-A leader 3-11 Tetramer-VMAPRTLVL-BV421	50 tests	¥209,000