

販売中止のご案内

**BioGenex 社商品 販売中止のお知らせ**

お客様各位

拝啓 時下ますますご清祥のこととお慶び申し上げます。

平素は当社商品に対しまして、格別のお引き立てを賜り、厚く御礼申し上げます。

このたび、BioGenex 社免疫組織化学染色 (IHC)、In Situ Hybridization (ISH) 関連商品につきまして、2020年3月をもちまして取り扱いを終了させていただきますので、謹んでご案内申し上げます。

お客様には大変ご迷惑をお掛けいたしますが、ご理解賜りますようお願い申し上げます。

敬具

2019年10月

販売中止日：2020年3月31日

対象品：

Code No.	商品名	クローン名	包装
B-AM4875M	Anti-ACTH mAb	AH26	6 mL
B-AM0905M	Anti-Actin, Muscle Specific mAb	HHF35	6 mL
B-AM1285M	Anti-Actin, Smooth Muscle mAb	1A4	6 mL
B-MU256UC	Anti-Androgen Receptor mAb	F39.4.1	1 mL
B-PU347UP	Anti-Bax Protein pAb	Polyclonal	1 mL
B-AM1355M	Anti-BCA225 mAb	CU18	6 mL
B-AN7235M	Anti-Bcl-2 mAb	EP36	6 mL
B-AM2875M	Anti-Bcl-2 Oncoprotein mAb	bcl-2/100	6 mL
B-AM7085M	Anti-Bcl-6 mAb	LN22	6 mL
B-AN7785M	Anti-Beta-Catenin mAb	EP35	6 mL
B-NU778UC	Anti-Beta-Catenin mAb	EP35	1 mL
B-MU178UC	Anti-Beta-Tubulin IV mAb	ONS1A6	1 mL
B-MU122UC	Anti-Beta-Tubulin mAb	DM-1B	1 mL
B-AM3595M	Anti-bFGF mAb	bFGF88	6 mL
B-AM2315M	Anti-BLA36 mAb	A27-42	6 mL
B-MU231UC	Anti-BLA36 mAb	A27-42	1 mL
B-MU429UC	Anti-CA125 mAb	0v185:1	1 mL
B-AM3325M	Anti-Caldesmon, High MW, Smooth Muscle mAb	h-CD	6 mL
B-AM3335M	Anti-Calponin mAb	CALP	6 mL
B-MU333UC	Anti-Calponin mAb	CALP	1 mL
B-AN8215M	Anti-Calponin-1 mAb	EP63	6 mL
B-NU821UC	Anti-Calponin-1 mAb	EP63	1 mL
B-AM5835M	Anti-Calretinin mAb	2E7	6 mL
B-PU413UP	Anti-Calretinin pAb	Polyclonal	1 mL
B-AM4905M	Anti-CD1a mAb	010	6 mL
B-MU490UC	Anti-CD1a mAb	010	1 mL

製造販売元

**MBL** 株式会社 医学生物学研究所

A JSR Life Sciences Company <https://ivd.mbl.co.jp/>

<お問い合わせ先>

〒105-0012

東京都港区芝大門2丁目11番8号 住友不動産 芝大門二丁目ビル

TEL : 03-6854-3611 (代表) FAX : 03-6854-3615

E-mail : kensa@mbl.co.jp

Code No.	商品名	クローン名	包装
B-AM4215M	Anti-CD4 mAb	4B12	6 mL
B-MU421UC	Anti-CD4 mAb	4B12	1 mL
B-AN7225M	Anti-CD4 mAb	EP204	6 mL
B-AM4305M	Anti-CD5 mAb	4C7	6 mL
B-MU430UC	Anti-CD5 mAb	4C7	1 mL
B-AM7025M	Anti-CD7 mAb	LP15	6 mL
B-AM2615M	Anti-CD8 mAb	T8	6 mL
B-AN7405M	Anti-CD8 mAb	SP16	6 mL
B-AN8515M	Anti-CD11b/ITAM mAb	EP45	6 mL
B-AN8225M	Anti-CD11c mAb	EP157	6 mL
B-NU822UC	Anti-CD11c mAb	EP157	1 mL
B-AN8145M	Anti-CD14 mAb	EP128	6 mL
B-NU814UC	Anti-CD14 mAb	EP128	1 mL
B-MU302UC	Anti-CD15 mAb	BRA4F1	1 mL
B-AN7295M	Anti-CD19 mAb	EP169	6 mL
B-MU439UC	Anti-CD22 mAb	FPC1	1 mL
B-AM2985M	Anti-CD29 mAb	JB1a	6 mL
B-AM2325M	Anti-CD31 (Endothelial Cell) mAb	JC/70A	6 mL
B-AM2415M	Anti-CD31 (PECAM-1) mAb	9G11	6 mL
B-MU241UC	Anti-CD31 (PECAM-1) mAb	9G11	1 mL
B-AM2365M	Anti-CD34 (Endothelial Cell) mAb	QBend/10	6 mL
B-AM4315M	Anti-CD35 mAb	RLB25	6 mL
B-NU769UC	Anti-CD38 mAb	SP149	1 mL
B-AN7325M	Anti-CD41 mAb	EP178	6 mL
B-NU732UC	Anti-CD41 mAb	EP178	1 mL
B-AM3055M	Anti-CD43 (T cell, Leukosialin) mAb	DFT-1	6 mL
B-AM3105M	Anti-CD44 mAb	DF1485	6 mL
B-AM2685M	Anti-CD56 (Natural Killer Cell, NCAM) mAb	NKH-1	6 mL
B-AM3145M	Anti-CD57 (Natural Killer Cell) mAb	NK-1	6 mL
B-AM3255M	Anti-CD66 mAb	BY114	6 mL
B-MU416UC	Anti-CD68 mAb	KP1	1 mL
B-MU549UC	Anti-CD68 mAb	CD68/G2	1 mL
B-AM3545M	Anti-CD71 (Transferrin Receptor) mAb	H68.4	6 mL
B-MU354UC	Anti-CD71 (Transferrin Receptor) mAb	H68.4	1 mL
B-AN7575M	Anti-CD82 (KAI1) mAb	EP160	6 mL
B-NU757UC	Anti-CD82 (KAI1) mAb	EP160	1 mL
B-MU355UC	Anti-CD99 (MIC2 Antigen) mAb	HO36.1.1	1 mL
B-AR7595R	Anti-CD117 (c-Kit) pAb	Polyclonal	6 mL
B-AM4235M	Anti-CD117 mAb	T595	6 mL
B-AN8185M	Anti-CD117 mAb	EP10	6 mL
B-NU818UC	Anti-CD117 mAb	EP10	1 mL
B-AN8375M	Anti-CD138 mAb	EP201	6 mL
B-AN7165M	Anti-CD146 mAb	EP54	6 mL
B-AN7375M	Anti-CD205 mAb	EP176	6 mL
B-AM5345M	Anti-CD227 (MUCIN 1) mAb	VU-4H5	6 mL
B-AM3925M	Anti-CDX-2 mAb	CDX2-88	6 mL
B-MU392AU	Anti-CDX-2 mAb	CDX2-88	1 mL
B-AM1265M	Anti-Chromogranin A mAb	LK2H10	6 mL
B-AM3185M	Anti-c-myc Protein mAb	9E10	6 mL
B-MU167UC	Anti-Collagen III mAb	HWD1.1	1 mL
B-AM3795M	Anti-Collagen IV mAb	COL-94	6 mL
B-NU815UC	Anti-Cyclin D1 mAb	EP12	1 mL
B-AN8475M	Anti-Cytokeratin 5 mAb	EP24	6 mL
B-AM1425M	Anti-Cytokeratin 8 mAb	C51	6 mL
B-AM1325M	Anti-Cytokeratin 13 mAb	AE8	6 mL
B-AM1465M	Anti-Cytokeratin 14 mAb	LL002	6 mL
B-MU146UC	Anti-Cytokeratin 14 mAb	LL002	1 mL
B-MU133UC	Anti-Cytokeratin, High MW mAb	AE3	1 mL
B-AM2915M	Anti-Cytokeratin, High MW mAb	34 $\beta$ E12	6 mL

製造販売元

**MBL** 株式会社 医学生物学研究所

A JSR Life Sciences Company <https://ivd.mbl.co.jp/>

<お問い合わせ先>

〒105-0012

東京都港区芝大門2丁目11番8号 住友不動産 芝大門二丁目ビル

TEL : 03-6854-3611 (代表) FAX : 03-6854-3615

E-mail : kensa@mbl.co.jp

Code No.	商品名	クローン名	包装
B-AM2545M	Anti-Cytomegalovirus mAb	BM204	6 mL
B-MU254UC	Anti-Cytomegalovirus mAb	BM204	1 mL
B-AM5705M	Anti-DOG1 mAb	1. 1	6 mL
B-AM2445M	Anti-Dystrophin mAb	Dys2 (Dy8/6C5)	6 mL
B-AM2225M	Anti-EBV mAb	1108-1	6 mL
B-AM3905M	Anti-E-Cadherin mAb	36	6 mL
B-MU390UC	Anti-E-Cadherin mAb	36	1 mL
B-AN7255M	Anti-E-Cadherin mAb	EP6	6 mL
B-NU820UC	Anti-Ep-CAM mAb	EP155	1 mL
B-AN7825M	Anti-ERG (Est-Related Gene) mAb	EP111	6 mL
B-NU782UC	Anti-ERG (Est-Related Gene) mAb	EP111	1 mL
B-AM3375M	Anti-Factor XIII A mAb	E980.1	6 mL
B-MU337UC	Anti-Factor XIII A mAb	E980.1	1 mL
B-AR7985R	Anti-FLI1 pAb	Polyclonal	6 mL
B-AR7665R	Anti-FSH pAb	Polyclonal	6 mL
B-AN8565M	Anti-GCDFP-15 mAb	EP95	6 mL
B-NU856UC	Anti-GCDFP-15 mAb	EP95	1 mL
B-MU336UC	Anti-GLEPP-1mAb	5C11	1 mL
B-AM5055M	Anti-GLUT-1 mAb	SPM498	6 mL
B-MU505UC	Anti-GLUT-1 mAb	SPM498	1 mL
B-AM5395M	Anti-Glypican-3 mAb	GPC3-88	6 mL
B-MU539UC	Anti-Glypican-3 mAb	GPC3-88	1 mL
B-PU082UP	Anti-HBcAg pAb	Polyclonal	1 mL
B-AR0215R	Anti-Hemoglobin A pAb	Polyclonal	6 mL
B-AM3625M	Anti-HPV16 mAb	Cam Vir-1	6 mL
B-MU374UC	Anti-J Chain mAb	JC88	1 mL
B-MU410UC	Anti-Ki67 mAb	K-2	1 mL
B-AN7875M	Anti-Luteinizing Hormone (LH) mAb	SP132	6 mL
B-AN8345M	Anti-MCM2 mAb	EP40	6 mL
B-MU361UC	Anti-Melan-A (MART-1) mAb	A103	1 mL
B-AM0015M	Anti-Melanoma mAb	HMB45	6 mL
B-MU001AU	Anti-Melanoma mAb	HMB45	1 mL
B-AM0775M	Anti-Melanoma-Associated Antigen mAb	NKI/C3	6 mL
B-MU077UC	Anti-Melanoma-Associated Antigen mAb	NKI/C3	1 mL
B-AM4335M	Anti-Mesothelin mAb	5B2	6 mL
B-MU213UC	Anti-Mitochondrial Antigen mAb	113-1	1 mL
B-AM7035M	Anti-MLH-1 mAb	ES05	6 mL
B-MU703UC	Anti-MLH-1 mAb	ES05	1 mL
B-NU816UC	Anti-MMP9 mAb	EP127	1 mL
B-AN7435M	Anti-MSH2 mAb	SP46	6 mL
B-NU743UC	Anti-MSH2 mAb	SP46	1 mL
B-AN8135M	Anti-MUC1 mAb	EP85	6 mL
B-NU813UC	Anti-MUC1 mAb	EP85	1 mL
B-AM4555M	Anti-MUC4 mAb	1G8	6 mL
B-MU455UC	Anti-MUC4 mAb	1G8	1 mL
B-AM4565M	Anti-MUC5AC mAb	45M1	6 mL
B-MU456UC	Anti-MUC5AC mAb	45M1	1 mL
B-AM3585M	Anti-Mucin-2 (MUC-2) mAb	CCP58	6 mL
B-MU358UC	Anti-Mucin-2 (MUC-2) mAb	CCP58	1 mL
B-NU750UC	Anti-Muml (IRF4) mAb	SP114	1 mL
B-AR4965R	Anti-Myeloperoxidase (MPO) pAb	Polyclonal	6 mL
B-AM4325M	Anti-Myf-4 mAb	L026	6 mL
B-MU432UC	Anti-Myf-4 mAb	L026	1 mL
B-AM0125M	Anti-Myoglobin mAb	MG-1	6 mL
B-MU109UC	Anti-Myosin, Skeletal Muscle mAb	MY-32	1 mL
B-AM7015M	Anti-Napsin A mAb	IP64	6 mL
B-AM0735M	Anti-Neurofilament mAb	NE-14	6 mL
B-AN7385M	Anti-NGFR mAb	EP31	6 mL
B-AN8305M	Anti-Oct 2 mAb	EP115	6 mL

製造販売元

**MBL** 株式会社 医学生物学研究所

A JSR Life Sciences Company <https://ivd.mbl.co.jp/>

<お問い合わせ先>

〒105-0012

東京都港区芝大門2丁目11番8号 住友不動産 芝大門二丁目ビル

TEL : 03-6854-3611 (代表) FAX : 03-6854-3615

E-mail : kensa@mbl.co.jp

Code No.	商品名	クローン名	包装
B-AN7245M	Anti-Oct 4 mAb	EP143	6 mL
B-AM5405M	Anti-p16 (INK4a) mAb	G175-405	6 mL
B-MU540UC	Anti-p16 (INK4a) mAb	G175-405	1 mL
B-AC6015M	Anti-p16 + Ki67 Cocktail mAb	G175-405 +	6 mL
B-MU434UC	Anti-p21/WAF1 mAb	4D10	1 mL
B-MU195UC	Anti-p53 Protein mAb	BP53-12-1	1 mL
B-AM2395M	Anti-p53 Protein mAb	D07	6 mL
B-MU239UC	Anti-p53 Protein mAb	D07	1 mL
B-AM2405M	Anti-p53 Protein mAb	1801	6 mL
B-MU240UC	Anti-p53 Protein mAb	1801	1 mL
B-AM4185M	Anti-p63 mAb	4A4	6 mL
B-MU418UC	Anti-p63 mAb	4A4	1 mL
B-AN4495M	Anti-P504S (AMACR) mAb	13H4	6 mL
B-NU449UC	Anti-P504S (AMACR) mAb	13H4	1 mL
B-AN5385M	Anti-P504S (AMACR) mAb	RBT-AMACR	6 mL
B-AM2525M	Anti-PCNA mAb	PC10	6 mL
B-PU376UP	Anti-PDGF (Platelet-Derived Growth Factor) pAb	Polyclonal	1 mL
B-AM3915M	Anti-P-Glycoprotein mAb	MDR88	6 mL
B-AR0405R	Anti-Placental Lactogen pAb	Polyclonal	6 mL
B-PU040UP	Anti-Placental Lactogen pAb	Polyclonal	1 mL
B-AM2285M	Anti-PLAP (Placental Alkaline Phosphatase) mAb	PL8-F6	6 mL
B-AN8445M	Anti-PMS2 mAb	EP51	6 mL
B-NU844UC	Anti-PMS2 mAb	EP51	1 mL
B-AM0135M	Anti-PSAP (Prostate Specific Acid Phosphatase) mAb	B01-94-21M-NA	6 mL
B-AN7145M	Anti-PSMA (Prostate Specific Membrane Antigen) mAb	EP192	6 mL
B-AN7465M	Anti-PTEN mAb	SP218	6 mL
B-AM5435M	Anti-RCC (Renal Cell Carcinoma) mAb	RCC-26	6 mL
B-AN8335M	Anti-SOX2 mAb	EP103	6 mL
B-AN8265M	Anti-Survivin mAb	EP119	6 mL
B-NU826UC	Anti-Survivin mAb	EP119	1 mL
B-AM3635M	Anti-Synaptophysin mAb	Snp88	6 mL
B-MU363UC	Anti-Synaptophysin mAb	Snp88	1 mL
B-AM4595M	Anti-Tau mAb	Tau-5	6 mL
B-MU377UC	Anti-TGF alpha mAb	TGF88	1 mL
B-AM034A5	Anti-Thyroxine mAb	D5	6 mL
B-AM5295M	Anti-TIA-1 mAb	2G9A10F5	6 mL
B-MU529UC	Anti-TIA-1 mAb	2G9A10F5	1 mL
B-AR1255R	Anti-Toxoplasma gondii pAb	Polyclonal	6 mL
B-PU125UP	Anti-Toxoplasma gondii pAb	Polyclonal	1 mL
B-AM0335M	Anti-TSH (Thyroid Stimulating Hormone) mAb	5404	6 mL
B-AR4835R	Anti-VEGF pAb	Polyclonal	6 mL
B-PU483UP	Anti-VEGF pAb	Polyclonal	1 mL
B-PU530UP	Anti-VIP pAb	Polyclonal	1 mL
B-HK0575K	Antigen Retrieval AR-10 Solution (10X Concentrated)	-	100 mL
B-HK05820	Antigen Retrieval AR-10 Solution (Ready-to-Use)	-	1000 mL
B-HK0585K	Antigen Retrieval AR-10 Solution (Ready-to-Use)	-	250 mL
B-HK0805K	Antigen Retrieval Citra Plus Solution (10X Concentrated)	-	100 mL
B-HK0809K	Antigen Retrieval Citra Plus Solution (10X Concentrated)	-	500 mL
B-HK08120	Antigen Retrieval Citra Plus Solution (Ready-to-Use)	-	1000 mL
B-HK0815K	Antigen Retrieval Citra Plus Solution (Ready-to-Use)	-	250 mL
B-HK0865K	Antigen Retrieval Citra Solution (10X Concentration), 100 mL	-	100 mL
B-HK0869K	Antigen Retrieval Citra Solution (10X Concentration), 500 mL	-	500 mL
B-HK0875K	Antigen Retrieval Citra Solution (ready-to-use)	-	250 mL
B-HK08720	Antigen Retrieval Citra Solution (ready-to-use), 1000 mL	-	1000 mL
B-HK1025K	Avidin/Biotin Blocking Kit (Ready-to-Use)	-	60 slides
B-XT134SL	Barrier Slides, 24 x 40 mm	-	72 slides
B-HK3405K	Biotinylated Anti-Immunoglobulins (Ready-to-Use)	-	6 mL
B-HK3409K	Biotinylated Anti-Immunoglobulins (Ready-to-Use)	-	50 mL
B-HK1565K	Concentrated-Reagent Diluents	-	100 mL

製造販売元

**MBL** 株式会社 医学生物学研究所

A JSR Life Sciences Company <https://ivd.mbl.co.jp/>

<お問い合わせ先>

〒105-0012

東京都港区芝大門2丁目11番8号 住友不動産 芝大門二丁目ビル

TEL : 03-6854-3611 (代表) FAX : 03-6854-3615

E-mail : kensa@mbl.co.jp

Code No.	商品名	クローン名	包装
B-HK94125	Enhanced Antibody Diluent (Ready-to-Use)	-	25 mL
B-HK941YA	Enhanced Antibody Diluent (Ready-to-Use)	-	100 mL
B-HK545X0	EZ-AR Common Solution (DeWax Solution) (5X Concentrated)	-	1000 mL
B-HK546XA	EZ-AR Elegance 1 Solution (Citra-based) (Ready-to-Use)	-	1000 mL
B-HK547XA	EZ-AR Elegance 2 Solution (EDTA-based) (Ready-to-Use)	-	1000 mL
B-HK5845K	EZ-DeWax Solution, Concentrated	-	500 mL
B-HK1009K	Hematoxylin (Mayer)	-	250 mL
B-HK8735K	Nucleic Acid Retrieval Solution 1 (NAR-1)	-	250 mL
B-HK3305K	Peroxidase-Conjugated Streptavidin (Ready-to-Use)	-	6 mL
B-HK3309K	Peroxidase-Conjugated Streptavidin (Ready-to-Use)	-	50 mL
B-HK1115K	Peroxide Block	-	6 mL
B-HK11150	Peroxide Block	-	50 mL
B-HK0855K	Power Block (10X Concentrated)	-	100 mL
B-EK0025K	Protease XXIV 3-Pack	-	150 tests
B-HK5835K	Super Sensitive Wash Buffer (20X Concentrated)	-	500 mL
B-HK542XA	Super Sensitive™ DAB (2 component DAB Pack)	-	1000 slides
B-QS20060	Super Sensitive™ Double Staining Polymer Detection Kit I	-	60 slides
B-QP9009L	Super Sensitive™ Multilink® HRP Detection Kit Large Volume	-	500 slides
B-QP300XA	Super Sensitive™ Multilink® HRP Detection Kit Mega Volume	-	1000 slides
B-QD610YA	Super Sensitive™ One-Step Polymer-HRP Detection System	-	200 slides
B-QD620YI	Super Sensitive™ One-Step Polymer-HRP IHC Detection System,	-	500 slides
B-DF40050	Super Sensitive™ One-Step Polymer-HRP ISH Detection System	-	50 slides
B-QD40060	Super Sensitive™ Polymer-HRP IHC Detection System/DAB	-	60 slides
B-QD420YI	Super Sensitive™ Polymer-HRP IHC Detection System/DAB Large	-	500 slides
B-QD430XA	Super Sensitive™ Polymer-HRP IHC Detection System/DAB Mega	-	1,000
B-HK0795K	SuperMount Permanent Mounting Medium	-	15 mL
B-HK0797K	SuperMount Permanent Mounting Medium	-	50 mL
B-EK0015K	Trypsin 3-Pack	-	150 tests
B-PR026	Alu II DNA Probe	-	25 slides
B-PR205	Epstein Barr Virus Encoded RNA Probe	-	25 slides
B-PR251	HPV14 Probe	-	25 slides
B-HM015A	Hsa-miR-15a Probe	-	25 slides
B-HM021	Hsa-miR-21 Probe	-	25 slides
B-HM031	Hsa-miR-31 Probe	-	25 slides
B-HM034A	Hsa-miR-34a Probe	-	25 slides
B-HM100	Hsa-miR-100 Probe	-	25 slides
B-HM126	Hsa-miR-126 Probe	-	25 slides
B-HM135B	Hsa-miR-135B Probe	-	25 slides
B-HM144	Hsa-miR-144 Probe	-	25 slides
B-HM155	Hsa-miR-155 Probe	-	25 slides
B-HM192	Hsa-miR-192 Probe	-	25 slides
B-HM195	Hsa-miR-195 Probe	-	25 slides
B-HM200A	Hsa-miR-200a Probe	-	25 slides
B-HM205	Hsa-miR-205 Probe	-	25 slides
B-HM214	Hsa-miR-214 Probe	-	25 slides
B-HM409-3	Hsa-miR-409-3P Probe	-	25 slides
B-PR214	Kappa Probe	-	25 slides
B-PR215	Lambda Probe	-	25 slides
B-PR032	miR-Scramble Probe	-	25 slides
B-PR031	miR-U6 Probe	-	25 slides

製造販売元

**MBL** 株式会社 医学生物学研究所

A JSR Life Sciences Company <https://ivd.mbl.co.jp/>

<お問い合わせ先>

〒105-0012

東京都港区芝大門2丁目11番8号 住友不動産 芝大門二丁目ビル

TEL : 03-6854-3611 (代表) FAX : 03-6854-3615

E-mail : kensa@mbl.co.jp